

Blood Pressure Trainer



ENGLISH

Directions for Use



Table of Contents

Laerdal Recommends.....	
Items Included.....	
Skills Taught.....	
Attachment of Blood Pressure Arm to SimPad.....	
Attachment of Blood Pressure Cuff to Arm.....	
Care and Maintenance	

The Blood Pressure Trainer is an accessory to the SimPad System that allows for the measurement of noninvasive blood pressure, brachial and radial pulse checks.

The manufacturing quality of this simulator should provide many sessions of training when reasonable care and maintenance are practiced.

Laerdal Recommends

For use with standard stethoscope.

Items Included:

- (1) Blood Pressure Trainer*
- (1) Carry Case
- (1) Blood Pressure Cuff

*Must be used with SimPad System (sold separately).

Skills Taught:

- The measurement of Noninvasive Blood Pressure, either auscultated or palpated
- Radial and Brachial Pulse Checks

Blood Pressure Trainer Connections –

Attachment of Blood Pressure Arm to SimPad–

Attach serial cable and clear pneumatic tubing to back of SimPad (Photo 1). If necessary use the adapter cable included with SimPad.



Photo 1

Attachment of Blood Pressure Cuff to Arm –

- a. Place supplied Blood Pressure Cuff on arm.
- b. Attach clear tubing on cuff to clear tube attached to arm. (Photo 2)




Photo 2

See SimPad Directions for Use for complete Blood Pressure operating procedures.

Laerdal Recommends

Do not insert needles into blood pressure arm.

See SimPad Directions for Use for additional connection and operating procedures.

Care and Maintenance –

1. Clean with mild soap and water; do not submerge trainer or parts in cleaning fluids or water. Apply a light coat of talcum powder to the face and chest skin to achieve skin-like feeling.
2. Use only on clean surface. Avoid felt tipped markers, ink pens, acetone, iodine or other staining products and avoid placing the manikin on newsprint or inked lines of any kind.
3. To ensure longevity, trainer should be cleaned after each training session and a general inspection should be conducted regularly.
4. Modules and all other parts should be drained and air-dried thoroughly before storage and disinfected when needed.
5. Articulating parts will benefit from a light application of talcum powder prior to training sessions.
6. Store properly between teaching sessions.