

Optional Equipment/Accessories

Masks

85 15 00 Silicone Mask No. 00 85 16 00 Silicone Mask No. 0/1 85 17 00 Silicone Mask No.2 86 02 20 Silicone Child Mask 3-4 w/Multi-Function Mask Cover 87 02 20 Silicone Adult Mask 4-5+ w/Multi-Function Mask Cover 86 52 00 Multi-Function Mask cover 3-4 87 52 00 Multi-Function Mask cover 4-5

Patient Valves

56 02 00 Patient Valve 54 01 03 Lip Valve 54 01 05 Disk Membranes, pkg. 10 85 12 50 Patient Valve w/ Pressure Relief Valve 85 12 52 Pressure Relief Valve (35 cm H_2O)

Storage

For storage in small space, the Adult and Paediatric (but not the Preterm) ventilation bags may be folded as shown in fig. 1-2-3. Fold the Paediatric in same fashion as the Adult.

Long Term Storage

Laerdal Silicone Resuscitators and/or spare parts may be placed in long term storage. They should be periodically inspected and tested (at least yearly) according 3 to the Function Testing section in this manual.

Storage Pouch (Not illustrated)

The Adult, Paediatric or Preterm resuscitator can be stored dust proof in the transparent resealable pouch. A loop through a metal eyelet permits hanging the pouch.

Compact Case

The Adult, Paediatric or Preterm resuscitator and accessories can be stored as illustrated. The Compact Case can also be wall mounted.

Display Case (Not illustrated)

A fully assembled Adult, Paediatric or Preterm resuscitator can be stored for immediate use. The Display Case may be hung on optional Wall Mount.

Wall Bracket

Separate Wall Bracket is available

Hanging Loop

The resuscitator can be hung ready for immediate use. For mounting Hanging Loop on the resuscitator, see separate instructions enclosed with the Hanging Loop.

85 13 50 Patient Valve w/Press Relief Valve and Lock Clip 85 | | 03 Lock Clips, pkg. | 0

Ventilation Bags/Intake Valve

- 85 01 50 Preterm Bag, 240 ml 86 01 50 Paediatric Bag, 500 ml 87 01 50 Adult Bag, 1600 ml 51 01 12 O-rings, pkg. 10 87 54 00 Intake Valve 87 19 50 Flap Valves, pkg 3
- 51 04 04 Intake membranes, pkg. 10 51 01 03 Cap, pkg. 3

0

Adult

Reservoirs 53 19 01 O2 Reservoir 2.6 litres 55 19 01 O2 Reservoir, 0.6 litre

Containers

85 07 00 Display Case cpl., Preterm 86 03 00 Display Case cpl., Paediatric 87 06 00 Display Case cpl., Adult 86 04 10 Compact Case cpl. Preterm/Paediatric 86 04 20 Compact Case cpl. Adult

Technical Specifications

Performance data may vary a great deal with the conditions under which they were obtained. Consequently, the findings in one test are not directly comparable with data found in another test unless test conditions were identical.

The product is in compliance with the CE essential requirements of Council Directive 0434 93/42/EEC as amended by Council Directive 2007/47/EC

- Product meets the following Product Standards: EN/ISO 10651-4:2002, Lung ventilators-Particular requirements for operatorpowered resuscitators
- ISO 8382: 1988 Resuscitators intended for Use With Humans.
- ASTM F 920 93, Standard Specification for Minimum Performance and Safety Requierments for Resuscitators Intended for Use With Humans.
- AS 2488-1995 Resuscitators intended for Use With Humans.

Operating environmental limits:

Operating condition: -18;C to 60;C. (-0,4°F to 140°F) 15% to 95 %

Feasible oxygen concentration

The feasible O_2 concentrations are approximated values and depend on the O_2 concentration delivered.

ADULT:	Ventilation bag volume:	1600 m
	Reservoir bag volume:	2600 m

Delivered O₂ concentrations under various test conditions:

O ₂ -flow lpm	Tidal vol. (ml) × bag cycling rate per min.O ₂ -concentrations (%) using reservoir (without reservoir).					
	400×12	400×24	600×12	600×24	1000×12	1000x24
3	74 (38)	51 (39)	58 (34)	40 (34)	44 (33)	33 (30)
8	100 (44)	100 (44)	100 (40)	68 (40)	78 (38)	51 (34)

100 (51) 100 (50) 100 (47) 100 (47) 100 (42) 75 (36)



- 85 05 00 Expiration Diverter (OD 30 mm) 87 10 00 Silicone Extension Tube (28 cm) 85 09 00 Manometer Connector 87 04 00 Laerdal Head Strap w/Attachment Ring 87 13 00 Attachment Ring f/Standard Head Strap 87 01 20 Hanging Loop 51 17 00 Wall Bracket 52 11 00 Wall Mount, Paediatric/Preterm displ. case 57 20 00 Wall Mount, ad. displ. case 87 05 50 Wall Poster reassembly guide 87 09 50 Directions for Use 53 19 07 Intake Valve Outer Part (23mm OD)
- 53 04 00 Airways, set of 4

	Relative Humidity
Storage environmental	limits:
Storage :	-40¡C to 70¡C.
	(-40°F to 158°F)
	40% to 95 %
	Relative Humidity
Dead space of Patient	Valve

Dead space of Patient Valve: Approx. 7,0 ml for all models

Expiratory resistance: Approx. 2,6 cmH₂O Measured with airflow of 50 lpm

Inspiratory resistance: w/reservoir approx. 4,2 cmH₂O w/o reservoir approx. 3,1 cm H_2O Measured with airflow of 50 lpm

Attainable delivery volume						
Adult:	Approx. 800m					
Paediatric:	Approx. 320m					
Preterm:	Approx. 150m					

Test conditions: Compliance 0,02 l/cm H₂O, Resistance 20 cm H₂O/l/s No leakage; Pressure Relief Valve overridden.



PAEDIATRIC: Ventilation bag volume: 500 ml Reservoir bag volume: 600 ml

flow	Tidal vol. (ml) x bag cycling rate per min.O ₂ -concentrations (%) using reservoir (without reservoir).					
	20×40	20×60	150×20	150x30	300×12	300x24
	100 (97)	100 (97)	98 (56)	78 (57)	85 (48)	56 (46)
	100 (100)	100 (100)	100 (70)	100 (70)	100 (58)	100 (57)
	100 (100)	100 (100)	100 (82)	100 (83)	100 (71)	100 (70)

E	R	Μ	ŀ	

Ventilation bag volume: 240 ml Reservoir bag volume: 600 ml

flow	Tidal vol. (ml) x bag cycling rate per min.O ₂ -concentrations (%) using reservoir (without reservoir).					
	20x40	20×60				
	100 (98)	100 (97)				
	100 (100)	100 (100)				
	100 (100)	100 (100)				

Spontaneous breathing patient

low	Tidal vol. (ml) \times bag cycling rate per min.O ₂ -concentrations (%) using reservoir (without reservoir).					
	Adult	Paediatric Preterm				
	600×20	300×20	150×25	20×60	20×60	
	44 (39)	66 (49)	99 (62)	100 (99)	100 (99)	
	81 (54)	98 (62)	99 (75)	100 (100)	100 (100)	
	96 (74)	98 (79)	99 (87)	100 (100)	100 (100)	

Laerdal Silicone Resuscitator products, accessories and parts are carefully engineered and produced using materials that are suitable for the purpose. Care in accordance with these Directions for Use will help ensure that each product has a long and useful lifetime. (Tested in 100 cycle decontamination study)

Material Chart					
Article	Comp	oonent	Part	Material	
Storage Pouch				Polyethylene	PE
Compact Case	Case			Polypropylene	PP
	Partiti	on wall		Acryl nitrilbutadiene styrene	e ABS
Display Case	Case			Polypropylene	PP
	Wind	OW		Styrene acrylonitril	SAN
	Tray			Acryl nitrilbutadiene styrene	ABS
	Lock			Polyamide	PA
Bag	Bag			Silicone rubber	SI
Dug	Valve	Connector		Polysulfone	PSU
	O-Rin	σ		Fluorelastomer	VITON
Patient Valve	Unne	'δ r Housing		Polysulfone	PSU
	Patien	nt side Housing		Polysulfone	PSU
	Lin Va			Silicope rubber	SI
	Disk N	Membrane		Silicone rubber	SI
	Pross	ire Relief Valve	Stem	Polysulfone	PSLI
	110350		Housing	Polysulfone	PSU
			Spring	Staipless steel	150
			Buch	Silicopo nubbor	CI
Intako Valvo	Outor	. part	DUSIT	Polyculfono	DCI I
IIILANC VAIVE	Innor	part		Polysulfono	PCLI
	Cap	part		Polysulfone	PSU
	Cap Elap V	/ahua		Silicopo mubbon	r SO
	Fiap v	Manahana		Silicone rubber	51
	Intake	riembrane		Silicone rubber	21
O2 Reconvoir Bag	Pasan	voir Pag		Polyainyl chlorido	PV/C
OZ Neservoir bag	Coup	ling for bag		Polysulfone	PSLI
	Coup	ing for dag		Folysullone	130
Macks	No.00)_0/1_2		Silicone rubber	SI
1 IdSKS	No.00	A A 5+		Silicopo rubbor	51
	No.0-	10		Silicono rubbor	SI
	140.0-	1-2		Silicone rubbei	51
Optional Equipment					
Mask Cover				Polysulfone	PSU
Lock Clip				Stainless steel	
Head Strap w/Ring	Strap			Polyvinyl chloride	PVC
riedd od up wirting	Attack	hement Ring		Polycarbonate	PC
Expiration Diverter	Housi	ng		Polysulfone	PSU
	Cente	er øasket		Silicone rubber	SI
	Extern	nal gasket		Silicone rubber	SI
Extension Tube	Tube	iai gusitet		Silicone rubber	SI
	Coup	ling		Polysulfone	PSU
Manometer Connect	or			Polysulfone	PSU
Hanging Loop	.01			Silicope rubber	SI
Wall Mount				Acryl nitrilbutadiene styrene	ABS
Wall Bracket				Acetal	POM
				/ cetu	1011
Shipping weights a	nd dimensio	ons			
Cat. Nos.		Weights		Dimensio	ons
850050	340g	l2 oz		25 x 14,5 x 13 cm	$10 \times 5.7 \times 5.1$ in
850051	380g	13 oz		25 x 14,5 x 13 cm	$10 \times 5.7 \times 5.1$ in
850053	830g	lb 3 oz		24 x 15 x 16 cm	9.4 x 5.9 x 6.3 in
850055	1640g	3 lb 10 oz		37 x 33 x 12 cm	14.6 × 13 × 4.7 in
860050	370g	13 oz		25 x 14,5 x 13 cm	10 x 5.7 x 5.1 in
860051	510g	l lb 2 oz		25 x 14,5 x 13 cm	$10 \times 5.7 \times 5.1$ in
860052	440g	16 oz		25 x 14.5 x 13 cm	10 x 5.7 x 5.1 in
860053	920g	2 lb		24 x 15 x 16 cm	9.4 x 5.9 x 6.3 in
860055	1760g	3 lb 14 oz		37 x 33 x 12 cm	14.6 x 13 x 4.7 in
860056	3900	14 07		25 x 14.5 x 13 cm	$10 \times 5.7 \times 51$ in
870050	520g	lh 2 oz		$25 \times 145 \times 13$ cm	$10 \times 5.7 \times 51$ in
870051	700g	lh 9 07		$25 \times 145 \times 13$ cm	$10 \times 57 \times 51$ in
870052	625g	lb 6 07		$25 \times 145 \times 13$ cm	$10 \times 57 \times 51$ in
870052	1060a	7 lb 5 oz		$23 \times 16, 3 \times 15$ Cm	94 2 59 2 42 10
870055	2090g	4 lb 10 oz		27 x 13 x 10 ull 37 x 23 x 15 cm	$146 \times 12 \times 47$ in
010000	20708	10 10 OZ		21 X 22 X 12 CIII	1 T.U. X. 1 J. X. 4./ IN

The Laerdal oxygen reservoir bag used with this product contains DEHP.

The product is only intended for short-term/transient patient application and it is not directly, or through a liquid, in contact with the patient.

When cleaned and used according to these Directions for Use, there is no added risk involved when treating children or pregnant or nursing women.

ENGLISH Directions for Use

Laerdal Silicone Resuscitators



www.laerdal.com

© 2010 Laerdal Medical AS, All Rights Reserved. 4492 rev F Printed in Norway



Scope

This manual provides the information required to fully utilise the Laerdal" Silicone Resuscitator, and to help assure safe and trouble free operation over a maximum period of time. The Laerdal Silicone Resuscitator is available in three sizes: Adult, Paediatric and Preterm. The 3 sizes have many features and functions in common and are therefore described jointly whenever possible. The detailed description of design and function should be read carefully by every user of the Laerdal Silicone Resuscitators.

It is mandatory that anyone who uses a manual resuscitator receive adequate instruction. It is also helpful if the student practices bag-valve-mask ventilation on realistic training manikins, such as the Laerdal manikins.

Cautions and warnings

Read these Directions for Use carefully and become thoroughly familiar with the operation and maintenance of the Laerdal Silicone Resuscitator before using it.

- Resuscitators should only be used by persons who have received adequate training.
- Federal law (US) restricts this device to sale by or on the order of a physician.
- Resuscitators should not be used with supplemental oxygen where smoking is permitted or when fire, flame, oil or grease is in close proximity.
- Resuscitators should not be used in toxic or hazardous atmospheres.
- Before first time use of the resuscitator parts and its accessories, decontamination is necessary
- The use of third party products and oxygen delivery devices (e.g. filters and demand valves) with the Laerdal Silicone Resuscitator may have an affect on LSR performance. Please consult with the manufacturer of the third party device to verify compatability with the LSR and obtain information on possible LSR performance changes.
- Do not use parts other than genuine Laerdal parts. Use of non-Laerdal parts may affect safety and/or performance.
- Laerdal strongly discourages the use of rinsing and drying agents. Such agents may not be compatible with the materials used in the Laerdal Silicone Resuscitator and may affect the material and/or performance.

Service

Laerdal Silicone Resuscitators are designed and engineered for utility and economy. All components, parts and assemblies listed in the parts list may be replaced, when necessary, by the operator who is directed to carefully inspect them during decontamination procedures. Under normal conditions of use, no regularly scheduled factory service should be necessary. However, whenever a question arises we hope you will contact Laerdal Medical AS or an authorised Laerdal distributor.

Product specifications are subject to change without notice.

Limited warranty

Please refer to the Global Warranty statement for additional terms and conditions (www.laerdal. com)

Manufactured by:

LAERDAL MEDICAL AS P.O. Box 377, N-4002 Stavanger Tel. +47 51 51 17 00, Fax +47 51 52 35 57 E-mail: laerdal.norway@laerdal.no

Distributed in USA by: LAERDAL MEDICAL CORPORATION 167 Myers Corners Road, P.O. Box 1840 Wappingers Falls, New York 12590-8840 Tel. (800) 431-1055, +1 (845) 297-7770 Fax (800) 227-1143, +1 (845) 298-4545 E-mail: customerservice@laerdal.com

Distributed in Canada by: LAERDAL MEDICAL CANADA LTD. 151 Nashdene Rd., Unit #45 Toronto, ON, Canada, MIV 4C3 Tel. +1 (416) 298-9600 Fax +1 (416) 298-8016 E-mail: savelives@laerdal.ca



Indications for Use

The Laerdal Silicone Resuscitator is a self-inflating manual resuscitator that is intended for patients requiring total or intermittent ventilatory support. The Laerdal Silicone Resuscitator provides positive pressure ventilation and allows spontaneous breathing either with a 22 mm ID (inner diameter) facemask port, through an artificial airway or with a facemask that has a 15 mm OD (outer diameter) connection. The Preterm model is intended for patients below 2,5 kg (5,5 lb), the Paediatric model is intended for patients from 2,5 (5,5 lb) to 25 kg (55 lb), and the Adult model is intended for patients over 25 kg (55 lb).

Ventilation with ambient air

Resuscitator ventilation without supplemental oxygen is possible.

Ventilation with oxygen

The Laerdal Silicone Resuscitator can be connected to an O_2 source via the oxygen nipple. Concentrations delivered to the patient depend on O_2 flow rate, use (or non-use) of a Reservoir Bag and ventilation technique, e.g. tidal volume, ventilation frequency, time relations during compression-release cycles. See Technical Specifications.

Inhalation of supplemental oxygen

A patient who breathes spontaneously can inhale O_2 through the resuscitator with minimal resistance. Attachment of the reservoir increases O_2 concentration. See Technical Specifications.

The mask can be hand held or strapped to the face.

Safety when using oxygen

- I. Build-up and transfer of high pressure to the patient is prevented since excess O_2 is vented to atmosphere over the outlet membrane of the Intake Valve.
- 2. When O_2 supply is insufficient, adequate ventilation volume is ensured by intake of ambient air over the intake membrane of the Intake Valve.
- 3. A Reservoir Bag that stays flat during the whole ventilation cycle is a visual indication that no, or little supplemental O_2 is being provided.

Pressure Relief Valve

The Preterm and Paediatric resuscitators feature a Patient Valve with a pressure limiting device mounted on the upper valve housing. If patient airway pressure exceeds 35 cm H_2O , the device opens to reduce

the risk of stomach distention and barotrauma. A hissing sound can be heard when the device opens.

When higher airway pressures are necessary, the operator can keep the Pressure Relief Valve closed with the tip of index finger while squeezing the bag.

A Lock Clip (optional) can be used as an alternative to finger pressure.

Accessories

Masks

The Laerdal Silicone Resuscitator can be combined with the following mask types and sizes:

a) Circular Infant Masks 00, 0/1, 2 b) Laerdal Child Mask 3-4 and Laerdal Adult Mask 4-5+

For difficult facial anatomies the Multi-Function Mask Cover is used to assist in getting a better mask seal.

All masks are transparent to enable the user to observe the patient's face and lip colour and the temporary fogging caused by exhalations.

Mask connection

The Patient Valve has a standard 15 (ID)/22(OD) mm patient port which connects to all standard masks or tube adapters. The Laerdal Masks 4-5+ and 3-4, plus the mask size 2, fit outside the patient valve connector. All other infant sizes fit inside, to reduce deadspace.

To use the Laerdal Head Strap











For Laerdal Child Mask 3-4 and Laerdal Adult Mask 4-5+, place the correct size Multi/Function Mask Cover over the mask connector. Fasten end of strap into the hooks on the cover. Tighten just enough to provide an airtight seal between mask and face.

For the Infant Mask 2, use the Attachment Ring supplied with the Laerdal Head Strap.

Expiration Diverter

An Expiration Diverter with two silicone gaskets can be snapped onto the Patient Valve.

The diverter provides an airtight seal to the valve housing but does not prevent the swivel function (possibility of horizontally rotating the bag without interfering with the position of mask or tube) of the valve connector. The diverter will provide an airtight seal when expired air is free flowing. The use of the Expiration Diverter with a restriction device (e.g. PEEP Valve) may cause some air leakage around the silicone gasket of the Expiration Diverter.

Equipment for measuring, scavenging or monitoring expired gases, can be attached to the standard (30 mm OD) outlet port of the diverter.

Manometer Connector

If used insert the Manometer Connector between the Patient Valve and the mask or tube adapter. Attach a manometer via tubing to the connector nipple (OD 6 mm) to monitor both inspiratory and expiratory pressures.

Extension Tube

The flexible Silicone Extension Tube may be used between ventilation bag and the Patient Valve. This extension tube makes it easier to ventilate when the patient is being transported. It also permits an operator to squeeze the bag against

a bed, strecher or themselves.



Practical Operation

a) When used in accordance with ISO 10651-4 the following resuscitator size recommendation applies: Adult for patients over 20 kg (44 lb), Paediatric for patients from 2.5 (5,5 lb) to 20 kg (44 lb) and Preterm for patients below 2,5 kg (5,5 lb).

When used to deliver tidal volumes as recommended by the AHA/ILCOR Guidelines 2000, the following applies. Adult for patients over 25 kg (5,5 lb), Paediatric for patients from 2,5 kg (5,5 lb) to 25 kg (55 lb) and Preterm for patients below 2,5 kg.

- b) Either connect the Patient Valve directly to the patient«s tube, or choose the appropriate size mask and attach it to the Patient Valve. Mask seal on difficult anatomies may be improved by using the Multi Function Mask Cover (Mask size 3-4 and 4-5+ only).
- c) Ventilate the patient by rhythmically compressing the bag for inspiration, allowing ample time between inspirations for patient's passive exhalation and bag re-expansion.
- d) Follow local guidelines for resuscitation.
- e) If the Patient Valve becomes contaminated with vomitus during ventilation, disconnect the resuscitator from the patient and clear the Patient Valve as follows.
- Tap the Patient Valve with the patient port against your gloved hand to shake free any contaminant and squeeze the silicone bag to deliver several sharp breaths through the Patient Valve to expel the contaminant.
- If contaminant does not clear; disassemble the Patient Valve and rinse.

Caution:

Visually inspect and test valve function to ensure proper operation of the Laerdal Silicone Resuscitator prior to patient use. Improper assembly of the flap valves, intake membrane, disk membrane and lip valve may affect performance. Misassembly of two lip valves may cause inadvertant EEP (End Expiratory Pressure) or prevent proper patient exhalation.

Decontamination

Thorough decontamination of the resuscitator components and accessories is necessary after each use. To reduce risk of cross contamination, follow steps below.

I. Washing and Rinsing

Washing and rinsing is always the first step in the decontamination process.

A Disassembly

- Disassemble the LSR into individual parts as shown in the Parts Illustration in Directions for Use, to make surfaces accessible to cleaning
- Separate Expiration Diverter (if used) into its three parts Separate Patient Valve into its four main parts
- For Preterm and Paediatric models, unscrew top of Pressure Relief Valve, but do not disassemble this part any further.
- Separate Intake Reservoir Valve into its six parts

IV. Sc

Inspection

4. Reassembly

I. Manual Cleaning B Rinse parts in a sink under cold running water from a tap. Submerge parts in warm tap water (30-40)/4 C / 86-104/4 F) ensuring that all surfaces are in contact with the warm water for at least 2 minutes before exposure to detergent. C Immerse all parts in hot tap water (60-701/4 C / 140-1581/4 F) containing a Dish Washing



CAUTION: Thorough cleaning and rinsing are the first and most important steps in the reprocessing of any reusable medical device. Without thorough cleaning and rinsing, it might not be possible to achieve high-level disinfection or sterilisation of the device.

I. The cycle has been validated on a Getinge Model A8666 validated to HTM2030. Ninhydrin protein detection test was used to qualify the process (to determine if any soil remained on the parts). Use of alternative washer/disinfector must be validated.

Meth

I. Ste (grav II. Ste (prev

> III. Cie (orth



CAUTION: Leave connectors in the necks of Ventilation Bags, Extension Tube, and Reservoir Bags during the entire decontamination procedure.

Laerdal strongly discourages the use of rinsing and drying agents. Such agents may not be compatible with the materials used in the Laerdal Silicone Resuscitator.

The use of non-validated cleaning, disinfecting or sterilisation methods may have adverse effects on the LSR material and/or performance.

Chose either the manual (I) or automatic (II & III) method below for cleaning the product.

- detergent ³.
- Thoroughly clean all surfaces using a brush as necessary. **D** Rinse all components free of detergent⁴ in warm tap water (30-401/4 C / 86-1041/4 F). Dry the components thoroughly ⁵
- **E** Inspect all components to confirm that they are CLEAN and DRY.
- II. Automatic Cleaning by Washer/Disinfector

Place parts in wire baskets. Cycle¹: 90-95¼C (194-203¼F) for more than 12 seconds. Total process time: approx. 52 min.² Use a Non-enzymatic alkaline detergent containing 2-5% NaOH3

III. Automatic Cleaning by Pasteurmatic Compact. ⁶ 30 min wash cycle at 32-43¹/₄C (90-110¹/₄F)

2. Includes pre rinse, main wash, rinse, final rinse and drying.

3. The Washing Detergent used in the validation - Olympic Chemicals Sprayclean 2000 (Non-enzymatic alkaline detergent containing 2-5% NaOH). Alternative detergents must be validated to show cleaning efficacy and material compatibility. Method has been validated using a common available tenside based Dish Washing Detergent (Zalo Ultra manufactured by Lilleborg AS). A pH neutral detergent solutions or hydrogen peroxide-based formulations may also be used for manual cleaning but must be validated to show they effectively clean the components. 4. CAUTION: If detergent or disinfectant residuals are allowed to dry on the resuscitator parts, surfaces may become

sticky, which may cause valve malfunction 5. Drying; the most effective method of drying is a fan assisted hot air cabinet, 50-70₁C (122-158¼F) for at least 30 minutes. Other drying methods may be used but must be validated to show they effectively dry the components. The Reservoir Bag must be dried by blowing air into the Reservoir Bag opening. 6. The cycle has been validated on a Pasteurmatic Compact from Olympic Medical.

2. Disinfection/Sterilisation

To obtain high-level disinfection/sterilisation of the resuscitator, the following 5 methods (I to V) have been validated and are recommended.

The sterilisation methods apply to all parts except reservoir bags, Head Straps, Wall Bracket, Storage Pouch and Containers. High-level disinfection methods apply to all parts.

Pasteurization applies to all parts except Wall Bracket, Storage Pouch and Containers

Method	Process paramet	Post-treatment				
	Parameters/Concentration	Exposure time				
terilisation						
I. Steam Autoclaving (gravity- displacement)	Autoclave at 132-137¼C (270 - 279°F)	15min. 00s (+ 30s)	Allow parts to cool and dry.			
II. Steam Autoclaving (prevacuum - pulse)	Autoclave at 134-137¼C (273-279¼F)	3min. 00s (+ 30s)				
ligh-level disinfection						
III. Cidex OPA (orthophtalaldehyde)	Conc.: 0,55% Ambient temperature	60 minutes	Remove traces of disinfectant by rinsing in warm tap water			
IV. Sodium Hypochlorite	Conc.: 0,5% Ambient temperature	20 minutes	(30-40¼ C / 86-104¼ F) for minimum 2 minutes. Dry the components thoroughly			
V. Pasteurization	Pasteurization cycle 70-75¼C (158-167¼F)	30 minutes	Dry the components thoroughly			

Carefully inspect all parts for signs of wear or damage. Worn or damaged components must be discarded and replaced with new components.

Reassemble resuscitator as shown in Parts/Assembly Illustration, in this Directions for Use.

Caution: Patient Valve reassembly

- Make sure that only one Lip Valve
- Cat.No. 54 01 03 is installed. If the valve housing does not tighten completely
- during reassembly, it may indicate that two lip valves
- have been mounted instead of one. Also, be sure not to mix the Disk Membrane for the Patient Valve
- with the Intake Membrane meant for the Intake Valve
- assembly. Test functions as described in Function Testing.

Intake Valve reassembly Reassembly as shown right.



Upper part of valve housing

- Valve housing, patient side

Disk membrane

Lip valve

Function Testing

Test valve functions to ensure proper operation of the resuscitator after each disassemblyreassembly. An O₂ Reservoir Bag is needed to complete the test procedures described below:

I. Intake/Reservoir Valve

- a) Compress the ventilation bag with one hand and close its neck opening with your other hand. Release the grip on the bag. Rapid bag reexpansion confirms efficient air intake.
- b) Close the neck opening and try to compress the bag. If the bag cannot be compressed with reasonable force, or if bag compression forces the air out between your hand and neck of the bag, the valve efficiently prevents backward leakage of air.

2.1 Patient Valve

a) Assure that a (single) Lip Valve has been installed in the Patient Valve. Attach the Patient Valve to the bag. Hold a Reservoir Bag over the patient port connector pressing with your thumb on the reservoir bag connector.

Ensure tight seal between the patient port and Reservoir Bag. Compress the bag with your other hand

several times. Inspect that the Lip Valve opens during compression.

Filling of the Reservoir Bag in this set-up confirms that the Patient Valve efficiently directs air to the patient.

b) With the filled Reservoir Bag held firmly to the valve connector, compress the Reservoir Bag while watching the external Disk Membrane.

Lifting of the Disk Membrane from its seat confirms that air is correctly directed to atmosphere instead of being returned to the ventilation bag.

2.2 Patient Valve with Pressure Relief Valve

Close patient port connector with your thumb while compressing the bag several times. Visual and audible opening of the relief valve confirms its operation.

3. Reservoir Flap Valves (located in the Intake Valve assembly.)

a) Do as described and shown in 2.1a above in order to fill the Reservoir Bag with ambient air. Attach reservoir to the Intake Valve and press on Reservoir Bag.

Compression of the Reservoir Bag and visual rise of the outlet Flap Valve confirms that the Reservoir Valve efficiently vents excessive gas to atmosphere.

b) Do as described and shown in 2.1a above in order to fill a Reservoir Bag with ambient air. Attach reservoir to the Intake Valve. With the Patient Valve in place and the reservoir attached to the Intake Valve, perform several compression-release cycles on the ventilation bag until the Reservoir Bag is flat and empty. Rapid reexpansion of the ventilation bag after flattening of the Reservoir Bag confirms that the Reservoir Valve efficiently lets in ambient air.













