
Thomas Tube Holder

ADULT

REF 600-10000 (qty.1)

REF 600-30000 (qty.1 - Green)

LATEX

LOT

Fits tube sizes 6.5mm (OD) - 21mm (OD) • Se adapta a tamaños de tubos de 6.5mm (OD) - 21mm (OD)
Para tubos de 6.5 mm (DE) a 21 mm (DE) • Per tubi da 6,5 mm (D.E.) -21 mm (D.E.) • S'adapte aux sondes de 6,5 mm à 21 mm (diamètre externe) • Passt zu den Schlauchgrößen 6,5 mm – 21 mm (jeweils äußerer Durchmesser) • Geschikt voor tubeafmetingen van 6,5 mm (buitendiameter) tot 21 mm (buitendiameter)
Passer til tubestørrelse 6,5 mm (OD) – 21 mm (OD) • Passar tubstorlek 6,5–21 mm (ytter diameter)
Passer til tubestørrelserne 6,5 mm (OD) – 21 mm (OD) • Sopii putkikokoihin 6,5–21 mm (ulkohalkaisija) •
• 合适的插管型号: 6.5mm (OD) - 21mm (OD) • Można stosować z rurkami intubacyjnymi o średnicy zewnętrznej od 6,5 mm do 21 mm

Non sterile • No es estéril • Não estéril • Non sterile • Non stérile • Unsteril • Niet steriel • Ikke steril •
Ikke-steril • Ikke-steril • Epästerili • 未消毒 • Produkt niesterylny

Thomas is a trademark of Laerdal Medical AS

www.laerdal.com

Laerdal
helping save lives

WARNING:

1. This device should only be used by trained personnel in accordance with the Directions for Use.
2. Single patient use.
3. Prior to use, inspect all parts to ensure they are free of damage and working properly.
4. Short term use only. Remove after 24 hours.
5. Patient should be constantly monitored while this product is in use.



CAUTION:

1. Federal law restricts this device to sale by or on the order of a physician or appropriately licensed practitioner.
2. Storage temperature: -34°C to +52°C (-30°F to 124°F)
3. Operating temperature: -18°C to +43°C (-0°F to 109.4°F)

USPatent # 5,402,776 # 5,513,633

Global warranty: See www.laerdal.com

Indications for use: Designed to secure single and double-lumen airway tubes and reduce the risk of accidental dislodgement.

1. Ensuring that the mouthpiece aperture is facing towards the patient's feet, slide the tube into V-wedge.
2. Advance bite block into patient's mouth ensuring lips are not caught between the teeth and bite block.
3. Run slide stick under patient's neck and then through the other side of holder.

4. Remove the slide stick and secure the fastening strap.
5. Tighten screw clamp securely against tube.

PLEASE NOTE:

If necessary, position removable hook tab on strap to accommodate smaller head or neck size.

It is important that the tube airway device is properly placed or positioned and the fastening strap secured before tightening the screw clamp.

PLEASE NOTE: Apply just enough pressure against the tube to hold securely in place. Use firm but sensible pressure. Overtightening may cause the tube to fully or partially close and prevent proper air flow.

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© 2015 Laerdal Medical AS. All rights reserved

Manufactured for:
LAERDAL MEDICAL AS
P.O. Box 377
Tanke Sivilandsgate 30
4002 Stavanger, Norway

Distributed in the USA by:
LAERDAL MEDICAL CORPORATION
167 Myers Corners Road, P.O. Box 1840
Wappingers Falls, New York 12590-8840
Tel. (800) 431-1055, +1 (845) 297-7770
Fax (800) 227-1143, +1 (845) 298-4545
E-mail: customerservice@laerdal.com

Distributed in Canada by:
LAERDAL MEDICAL CANADA LTD.
151 Nashdene Rd., Unit #45
Toronto, ON, Canada, M1V 4C3
Tel. (416) 298-9600
E-mail: savelives@laerdal.ca

Made in China

S181 rev J

www.laerdal.com