# Resusci Anne

Important Product Information

**Laerdal**  
helping save lives

English

## Resusci Anne QCPR / QCPR AED / QCPR D

---

### Notes

Important information about the product or its operation.

### Resusci Anne with Airway Head

The airways in the Airway Head cannot be completely sanitized, therefore, do not do:

- Mouth-to-mouth ventilation
- Mouth-to-mask ventilation
- Insertion of simulated vomit for suctioning.

### Resusci Anne AED

- When removing or replacing the chest skin, do not pull or damage the wires connecting the chest skin to the battery box.
- Use training pads with a white connector (Cat. No. 94 50 90).

### Care, Maintenance and Cleaning

For instructions on cleaning Resusci Anne faces, refer to the Cleaning Advice for Laerdal manikin faces insert. Keep the manikin clean and in a hygienic condition.

- Replace the airway assembly after every class, especially if mouth to mouth resuscitation is practiced.
- We recommend that you use a separate Laerdal Manikin Face for every student.
- If several students use one manikin face, thoroughly sanitize the manikin face after every use.
- Clean all skin parts regularly. Use warm soapy water or Laerdal manikin wipes.
- Use of Manikin Face Shields provides a clean barrier between your lips and the manikin face.

These can discolor the manikin:

- Pigments from lipstick and pens
- Latex gloves
- Using clothes other than Resusci Anne clothing.

### Warnings and Cautions

A Warning states a condition, hazard, or unsafe practice that can result in serious personal injury or death.

A Caution states a condition, hazard, or unsafe practice that can result in minor personal injury or damage to the manikin.

### Warnings

#### Resusci Anne QCPR D

- Observe all standard safety precautions for the use of defibrillators.
- Use a defibrillation chest skin when applying a defibrillator. Using another type of chest skin for defibrillation can damage the manikin.
- Do not defibrillate the manikin if the chest skin is not correctly fitted.
- To prevent overheating, only give a maximum of 2 x 360J defibrillator discharges per minute.
- Do not use cables or connectors with visible damage. Damaged cables and connectors can cause overheating, arcing and electrocution.
- Only connect the defibrillation paddles to the defibrillation connectors.
- Do not allow the manikin to contact electrically conductive surfaces or objects during defibrillation.
- Do not defibrillate the manikin in a flammable or oxygen-enriched atmosphere.
- Do not defibrillate the manikin when it is switched OFF or not functioning correctly.

### Cautions

Lubrication when using the Airway Head

- Use only lubricant provided by Laerdal Medical. Using silicon or other lubricant not approved by Laerdal can damage the airway.
- Lubricate instruments and tubes before insertion into the airway. It is difficult to insert non-lubricated instruments and tubes into the airway. Non-lubricated instruments and tubes can also damage the airway.

To prevent skin pitting on the Resusci Anne QCPR D manikin

- Do not apply conductive gel or conductive defibrillation pads intended for patient use.

### General

- Use of automatic chest compression machines may damage the manikin.

## Certification, Compliance and Labels

	The manikin is CE marked in accordance with Council Directive 1999/5/EC relating to Radio and Telecommunications Terminal Equipment (R&TTE), and Council Directive 2011/65/EU on Restriction on the use of certain hazardous substance (RoHS).
	Australian Radiocommunications and EMC Compliance Mark N25270
	CANADIAN ICES-003 STATEMENT CAN ICES-3 (B)/NMB-3(B)
	Manufacturer
	Date of Manufacture
	Dispose of in accordance with your country's recommendations.
	Reference order number
	Warning / Caution symbol

## FCC Statement

	This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: <ul style="list-style-type: none"><li>• This device may not cause harmful interference</li><li>• This device must accept any interference received, including interference that may cause undesired operation.</li></ul>
--	--

### Caution

Changes or modifications not expressly approved by Laerdal Medical could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Contains FCC ID: U3O-G2M5477

Contains FCC ID: VVXMLM400

## Specifications

Resusci Anne QCPR, QCPR AED, QCPR D  
Dimensions 69.7" x 20.5" x 9.8" (177 cm x 52 cm x 25 cm)  
Weight: 79.2 lbs (36 kg)  
Defibrillation Average of 2 x 360/j minute max (QCPR D only)  
Operation temperature: 32°F to 104°F (0°C to +40°C)  
Humidity: < 95% relative humidity  
Storage temperature: 5°F to +122°F (-15°C to +50°C)

## Ventilation

### Supported airway management tools

The following intubation devices have been tested on the Airway Head:

- LMA Classic 4
- LMA Classic 5
- LMA Unique 5
- LMA Fasttrack 4
- LMA Fasttrack 5
- Combitube 37 Fr
- LTS-D 4
- LTS-D 5
- Japanese Sumi NA

## Waste Handling

The manikin contains electronic components. Dispose of it at the applicable recycling facility in accordance with local regulations.

## European Directive 2012/19/EU (WEEE)

WEEE: this appliance is marked according to the European directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE). By ensuring this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste handling of this product.

The symbol on the product, or on the documents accompanying the product, indicates that this appliance may not be treated as household waste. Instead it shall be handed over to the applicable collection point for the recycling of electrical and electronic equipment. Disposal must be carried out in accordance with local environmental regulations for waste disposal. For more detailed information about treatment, recovery and recycling of this product, please contact your local city office, your household waste disposal service or the Laerdal representative where you purchased the product.

## Warranty

The Laerdal Resusci Anne QCPR / QCPR AED / QCPR D has a one-year limited Warranty. Refer to the Laerdal Global Warranty for terms and conditions.

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© 2013 Laerdal Medical AS. All rights reserved.