

# SimMan® Replaceable Parts



38 11 02  
Fastener for Neck Skin  
(pkg. 10)



200-03050  
Set, Pupils



38 11 07  
Teeth Upper Soft



38 11 06  
Teeth Upper Rigid  
(polyester)



38 11 05  
Neck Skin  
(pkg. 6)

211-10050  
Head Assembly

211-10150  
Head Skin

38 02 00  
Arm Assembly  
Blood Pressure

38 02 10  
Cuff Assembly  
Blood Pressure



300-00750  
Simulated  
Blood




38 07 00  
IV Arm, right



300-00650  
Manikin  
Lubricant

Stomach Bladder  
380484

31 20 29  
Skin & Veins, IV Arm



1005162  
Thigh Pads,  
set of 2

38 04 55  
Chest Skin  
with Chest Drain



211-11050  
Leg Assembly, right with  
Pedal Pulse

211-12050  
Leg Assembly, left with  
Pedal Pulse

1002961  
IV with Tubing  
IVR 5000



38 04 62  
Kit Adult Male Genitalia  
w/Valves

38 04 61  
Kit Adult Female Genitalia  
w/Valves





38 08 10  
Pad Set  
Shoulder and Hip



38 04 05  
Bladder Assembly Mid-Clavicular  
(Pneumothorax Chest)



38 04 06  
Bladder Set Mid-Auxillary  
(Pneumothorax Side)



200-00250  
Cricothyroid Membrane Tape  
(1 roll)



38 31 10  
Chest Tube Insertion Modules  
(pkg. 6)



38 04 10  
Post Set, ECG/defibrillation



245-16750  
Audio Cable

38 13 00  
Pulse Oximeter Probe



38 03 10  
Cable Assembly Manikin-Link Box  
(15 pin D-sub)




38 04 35  
Breathing Bladder



38 12 01  
Jet-Vent Adapter Kit



38 10 10  
Tubing Assembly Compressor-Manikin,  
Air/CO<sub>2</sub>



27 15 00  
Cable PC/Monitor to Link Box  
(9 pin D-sub)

### NOT ILLUSTRATED

- 1005146 Tubing Assy 40 ft  
Compressor-Manikin, Air/CO<sub>2</sub>
- 1005164 Cable Assy 25 ft  
Manikin-Link Box  
(15-pin D-sub)
- 260305 Power Cord (US) 110VAC
- 260306 Power Cord (EUR) 220VAC
- 260307 Power Cord (UK) 240VAC
- 380100 Link Box
- 380901 Jacket
- 380902 Trousers
- 211-19050 Directions for Use
- 211-09050 Backwards Compatibility Kit
- 380150 Patient Monitor SimMan
- 380201 Hardware Set, Arms to
- 380400 Chest Plate Assembly
- 380420 Lungs
- 380445 Shaft, Compression
- 380446 Foam in Chest
- 380460 Genitalia Blank
- 380470 Belly Plate Assembly
- 380471 Plug, Belly Plate
- 380475 Pin, Pelvis
- 380480 Thrust Assembly Cpl.

### Parts for authorized service personnel only

- 380300 Main Board
- 380350 Compression Board
- 380401 Speaker Set, Chest Plate
- 380402 Speaker Set, Mid-Auxillary
- 380411 Wire Harness, Chest Plate
- 380413 Cable Assy, Solenoid Valve
- 380430 Tubing Set w/Labels
- 380491 Pad Assy, Femoral Pulse Right
- 380492 Pad Assy, Femoral Pulse Left
- 380493 Cable Assy, Femoral Pulse
- 381000 Valve Assy, Airbank
- 381011 Connector, Female
- Dual-Distribution Panel Mount
- 381050 Manifold Assembly
- 381109 Upper Palate Assy
- 381110 Bladder Assembly, Head
- 381120 Shroud Assy w/Pulse-Sensors
- 381129 Foam Set, Carotid Pulse
- 381130 Bronchial Tree Cpl.
- 381133 Sock, Bronchial Tree
- 383100 Update Kit, Chest Drain



**Laerdal**  
helping save lives

[www.laerdal.com](http://www.laerdal.com)