



**Laerdal**<sup>®</sup>

helping save lives



## Preparing for Sudden Cardiac Arrest

### AED Trainer 3

A fast, accurate response from medical personnel when treating the cardiac arrest patient is essential if successful patient outcomes are to be achieved.

The AED Trainer 3 reflects the unique functionality of the Philips FR3 AED and provides a realistic training experience for medical staff to deliver optimal therapy when facing cardiac arrest.



# AED Trainer 3

## PHYSICAL

Weight 1.3 lbs (600g) (with batteries installed)  
Size HxWxD 8.6" x 5.2" x 2.2"  
(218 mm x 133 mm x 57 mm)

## OPERATION

Operating System On Data Card - 1 GB SD

Data Card also contains configuration software (multi-language) required to change device parameters.

All localization and software updates available from Laerdal download site.

## BATTERY

### AED Trainer 3

Type 4 AA Alkaline (any brand)  
Capacity >10 hours

### Remote Control

Battery Type Lithium CR2025 3V

## TRAINING PADS

Type HeartStart Training Pads III or Laerdal LINK Training Pads

## MATERIALS

Housing Acrylonitrile Butadiene Styrene (ABS)  
Control buttons Silicone  
Foot Polyurethane (PU)  
Infant/Child key Polycarbonate and Polybutylene Terephthalate (PC + PBT)  
Tether Silicone

## EMC CLASSIFICATION

The product is in compliance with the essential requirements of Directive 2004/108/EC on electromagnetic compatibility through compliance with:

- IEC 61000-6-3 Electromagnetic Compatibility (EMC) - Part 6-3:
- IEC 61000-6-1 Electromagnetic Compatibility (EMC) - Part 6-1:

(See DFU for additional electromagnetic conformity information)

## REGULATORY APPROVALS

CE mark, FCC Part 15 Class B Statement (US), CSA Certification (Canada), C Tick (N25270) (Australia)

## ENVIRONMENTAL CONDITIONS

Storage Temperature 0° - 40° C (32° - 104° F)  
Operating Temperature 10° - 35° C (50° - 95° F)  
Relative Humidity 0 - 90% non-condensing

## DISPOSAL

AED Trainer 3 contains electronic components which should be disposed of in accordance with local regulations.

## WARRANTY

One-Year limited. Please see Laerdal Global Warranty for terms and conditions.

## SERVICE

Repairs carried out by Laerdal Technical Centers.

## ORDERING INFORMATION

Part No.	Description
198-00150	AED Trainer 3, DFU, Quick Start Guide, 4 AA Alkaline Batteries, Training Infant/Child Key with Tether, Carry Case, Operating System Card, HeartStart Training Pads III, HeartStart Interconnect Cable, and External Manikin Adapter.

## Consumables

Part No.	Description
989803150181	HeartStart Training Pads III
989803150201	HeartStart Training Pads III Interconnect Cable
945090	Laerdal LINK Training Pads

## Accessories

Part No.	Description
198-00350	Remote Control (with CR2025 3V battery)
198-10050	Hard Carry Case
198-10250	1 GB Operating System Card w/Software
M5089A	External Manikin Adapters (5 pack)

Visit [www.laerdal.com/downloads](http://www.laerdal.com/downloads) for multi-language DFU, latest localizations and configuration software



Shown with Remote Control, not included with AED Trainer 3

877-LAERDAL (523-7325)  
[www.laerdal.com](http://www.laerdal.com)



**Laerdal**<sup>®</sup>  
helping save lives