


Parts Catalogue



Laerdal Parts Catalogue

October 2008 - Volume 4



Laerdal
helping save lives

Parts

The following pages are intended as a guide for identifying parts you may need to replace. Please note that the breakdown of a product does not specify the complete product (manikin, defibrillator, airway device, suction unit), but shows available spare parts and accessories.

If you need further information regarding the parts for your Laerdal products, please contact us. Trained personnel are available to answer your questions and provide additional technical data.

Spare parts for defibrillators (HeartStart) are not available for sale in the US through Laerdal Medical.

All parts listed for a specific product might not be available in all markets

*Product specifications are subject to change without notice.
Actual colours on the parts might differ from the colours in this catalogue.*

For online purchase please visit
www.laerdal.com



Training

Resusci Anne Modular System - Basic / SkillGuide	4-5
Resusci Anne Modular System - SkillReporter	6-7
Resusci Anne Modular System - CPR-D.....	8-9
Resusci Anne Modular System - AED SkillGuide.....	10-11
Resusci Anne Modular System - AED SkillReporter....	12-13
Resusci Anne Modular System.....	14-15
Resusci Junior	16-17
Resusci Baby	18-19
Little Anne - AED Little Anne.....	20
Little Junior	21
Baby Anne.....	22
ALS SkillMaster 4000	24-25
ALS SkillTrainer 200 (discontinued model)	26-27
Laerdal ALS Baby Trainer.....	28
Laerdal Infant Airway Management Trainer.....	29
Laerdal Airway Management Trainer	30
Laerdal IV Trainers.....	31
Resusci Anne Advanced SkillTrainer	32-33
Resusci Anne Simulator	34-35
SimMan.....	36-37
SimMan (discontinued model)	38-39
SimBaby	40-41
SimNewB	42-43
ALS Simulator	44-45
MegaCode Kelly Basic	46-47
MegaCode Kelly Advanced	48-49
ECG Kid	50
MegaCode Kid	51
Nursing Kelly	52-55
Nursing Anne	56-59
VitalSim	59
Nursing Kid	60-61
Nursing Baby	62-63
Neonatal Resuscitation Baby.....	64
Pneumothorax Trainer.....	65
Convalescent Kelly	66
NG Tube and Trach Care Trainer.....	67
Cath. and Enema Trainer	68
Virtual IV / Infant Virtual IV	69
Ultimate Hurt.....	70-71
First Aid Training Manikin	72
Tuff Kelly	73
Extri Kelly.....	74
AT Kelly Torso	75


Therapy

HeartStart MRx	78-81
HeartStart FR – FR2+	82
HeartStart HSI-FRx.....	83-84
HeartStart 4000.....	85
HeartStart FR (discontinued model).....	86
HeartStart 911 (discontinued model).....	87
HeartStart 3000/3000QR (discontinued models).....	88-89
Laerdal Silicone Resuscitators	90-91
Laerdal Suction Unit.....	92-93
Laerdal Compact Suction Unit.....	94-95
Laerdal Premier Suction Unit (discontinued model).....	96
Laerdal Suction Unit (discontinued model).....	97

Resusci Anne Modular System Cat. Nos. 31000x, 31001x, 31003x, 31004x



310200 Head cpl. basic



NOT ILLUSTRATED


042000 Carrying case Torso
183910 CPR training mat
311500 Directions for use
311590 Technical manual
152007 Paper bags f/face skin (pkg. 1000)

Parts for authorized service personnel are listed in the technical manual


153800 Side connection unit (SkillGuide manikin)




153900 SkillGuide cpl.




Use 4 x 1.5V size C batteries



Resusci Anne Modular System Cat. Nos. 31002x, 31005x



310400 Head cpl. w/electronics




NOT ILLUSTRATED


042000	Carrying case Torso
183910	CPR training mat
311500	Directions for use
311590	Technical manual
152007	Paper bags f/face skin (pkg. 1000)

Parts for authorized service personnel are listed in the technical manual


315000xx SkillReporter cpl.



Resusci Anne Modular System Cat. Nos. 32002x, 32005x



320100 Head cpl. w/automatic pulse



NOT ILLUSTRATED

042000	Carrying case Torso
183910	CPR training mat
311500	Directions for use
311590	Technical manual
152007	Paper bags f/face skin (pkg. 1000)
320256	Power resistor 150R 50W
291207	Cable main board - loadbox

Parts for authorized service personnel are listed in the technical manual

945004
Flat cable to manikin
(For MRx,
HS 4000, FR and FR2)



315000xx SkillReporter cpl.



325000 MicroHeartsim cpl.



321040
Barrel to manikin
(For MRx, HS 4000)



310681
Strap f/battery holder



310680
Battery holder
(Use 8 x 1.5V
size D batteries)



310650
Manikin keyboard



32 50 05
Cable main board




310500
Sensor unit cpl.



Resusci Anne Modular System Cat. Nos. 320060, 320080



310200 Head cpl. basic




NOT ILLUSTRATED

- 042000 Carrying case Torso
- 183910 CPR training mat
- 311500 Directions for use
- 311590 Technical manual
- 152007 Paper bags f/face skin (pkg. 1000)

Parts for authorized service personnel are listed in the technical manual

153900 SkillGuide cpl.



153800 Side connection unit (SkillGuide manikin)



110605 Blind plate



153620 Cable f/ hand position circuit board

153610 Circuit board f/ hand position unit

110604 Cover f/hand position unit



945050 Remote control for AED Trainer 2



945020 Training battery



945030 Softpack



945090 Training pads, Laerdal Link Technology


NOT ILLUSTRATED

- 945021 Label status AED Trainer 2 (pkg. 10)
- 945025 Board AED Trainer 2
- 945026 Speaker with cable
- 945027 Cable AED Trainer 2
- 945055 Programming kit for AED Trainer 2 (optional software and cable)
- 945060 Technical manual

Resusci Anne Modular System Cat. Nos. 320070, 320090



310400 Head cpl. w/electronics



NOT ILLUSTRATED

042000 Carrying case Torso
183910 CPR training mat
311500 Directions for use
311590 Technical manual
152007 Paper bags f/face skin (pkg. 1000)

Parts for authorized service personnel are listed in the technical manual

315000xx SkillReporter cpl.



310680 Battery holder (Use 8 x 1.5 V size D batteries)



NOT ILLUSTRATED

945021 Label status
945025 AED Trainer 2 (pkg. 10)
945025 Board AED Trainer 2
945026 Speaker with cable
945027 Cable AED Trainer 2
945055 Programming kit for AED Trainer 2 (optional software and cable)
945060 Technical manual

Resusci Anne Modular System

The Resusci Anne Modular System allows you to tailor make your manikin to meet the diversity of training needs, from the simplest CPR Torso to the full body / full featured CPR-D and First Aid versions.

Each Torso model can be equipped with various types of arms and legs.



Arms with soft legs



310300

Rescue module



312000

First Aid/Trauma module



312050

Bleeding control leg



090050

IV arm




312025

Simple arm and leg attachments allow you to build the manikin of your choice within seconds. No tools are required.

Update kits Resusci Anne Modular System

Update your Resusci Anne Basic
to SkillGuide Resusci Anne

154000 SkillGuide update kit




Update your Resusci Anne Basic
or SkillGuide to SkillReporter
Resusci Anne

311000 SkillReporter update kit




Update your Resusci Anne Basic or SkillGuide to CPR-D Resusci Anne

321000	Defib. update kit Basic/SkillGuide
--------	---------------------------------------





Update your SkillReporter Resusci Anne
to CPR-D Resusci Anne

321010	Defib. update kit SkillReporter
--------	------------------------------------






Additional parts for 180011 only




Parts for discontinued models

180001 / 180002 / 180003 only:




Parts for Resusci Junior models

180010 / 180011 / 180012 before G2005 update:



I43800 Head, complete



I44500 Update kit Basic to SkillGuide



NOT ILLUSTRATED


- I44900 Directions for use Resusci Baby
- I44910 Technical manual

Parts for authorized service personnel are listed in the technical manual

OLDER STYLE

- I40700 Chest skin

Additional parts for I4001 I only



I44512
Chest plate w/sen-
sor



I44511
Side connection
board for
foil-switch

I43907
Circuit board,
vent./comp.



I44550
Chest plate basic

I44490 Update Kit Guidelines 2005 (enable deeper compressions)









NOT ILLUSTRATED

050600	Soft pack
050610	Bag 4-pack
050710	Rib plate (pkg. 10)
020780	Family pack

For online purchase please visit
www.laerdal.com





291240
Valve board cpl.

290603
Airway valve base

290602
Airway valve housing



263020
Heart Rhythm CD Rom



291230
UMR cpl. w/electronics

291233
Internal cables

945004
Flat cable to manikin
(For MRx,
HS 4000, FR and FR2)



321040
Barrel to manikin
(For MRx, HS 4000)



110105
Wheeling
handle



292200 Carrying case



110104
Wheels, set of 2

290230 Chest skin with
4-connectors ECG & loadbox SkillMaster



290210
Chest skin

260030 Heartsim 4000 cpl.



263005
Link box



263010
CD Rom



263001
Remote control
(Use 4 x 1.5 V
size AAA batter-
ies)



271500
9-pin computer
cable



260306
Power cord (EU)

260305
Power cord (US)

260307
Power cord (UK)

NOT ILLUSTRATED

291830 Directions for use
291835 Technical manual
290215 Power resistor 50R 50W
293100 Cleaning kit

Parts for authorized service personnel
are listed in the technical manual

ALS SkillTrainer 200 Cat. Nos. 290020, 290021




945004
Flat cable to manikin
(For MRx,
HS 4000, FR and FR2)



321040
Barrel to manikin
(For MRx, HS 4000)



290203 Chest skin cpl.



290210
Chest skin ALS, no
electronics

260010 Heartsim 200



261000
Soft pack for
Heartsim 200

261020
Battery case lid

Use 6 x 1.5 V size AA batteries

NOT ILLUSTRATED

291805	Directions for use
291835	Technical manual
290215	Power resistor 50R 50W
293100	Cleaning kit

110105
Wheeling
handle



292200 Carrying case



110104
Wheels, set of 2

Laerdal ALS Baby Trainer Cat. Nos. 080030, 080031, 080035

083300 Head cpl.



083105 Pulse bulb




260010 Heartsim 200



Use 6 x 1.5V size AA batteries

141600 Compression slider
141700 Compression spring assembly



082010 Chest plate



082022 Stomach valves (pkg. 10)



083450 Chest skin cpl. w/ECG G2005




083050 ALS Baby Update kit Guidelines 2005
(enable deeper compressions)




142600 Pants



For 080035 only



080015 Intraosseous Trainer



082305 Leg replacement pads (pkg. 5)



144002 Tray for case



144010 Carrying case





NOT ILLUSTRATED


081500 Cleaning kit
081650 Directions for use
082025 Lungs/stomach (pkg. 6)

OLDER STYLE

083400 Chest skin cpl. w/ECG
083420 Chest skin (plain)



Laerdal Airway Management Trainer Cat. No. 250000



NOT ILLUSTRATED

250100	Foam block holder
253000	Directions for use
252600	Cleaning kit
250200	Head protector

090019 IV Torso, Basic

092010
IV Torso
inner shell

092003
Femoral
replacement pad

092012
IV Torso pulse tube
set

092001
Neck replacement
pad

092011
IV Torso
outer shell

110100
Carrying case

090022 IV Foot


Discontinued model

092210
IV Foot outer shell

092201
Replacement
pads for foot
(pkg. 5)

092211
IV Foot base

151-11050
Head cpl.



150-11150
Pupils (pkg.6)



310321
Fastening bolt for arm
(pkg.2)



151-10150
Arm right plain female



150-10150
IV Arm left female



201500
Compression spring



325-00350 Kit, Skin-
Vein Adult Female
Multi-Venous IV Arm



310310
Lower body



310345
Soft stomach



151-14050 Chest skin complete



151-12001
Remote Control FCC-ver.
151-12033
Remote Control CE-ver.



153910
Protecting pouch





150-19050
Software CD w/USB cable



310680
Battery holder
(Use 8x 1.5V size D batteries)

NOT ILLUSTRATED

310350	Clothing
150-11250	Upper airway cpl. wo/teeth
150-11450	Teeth set
150-11550	Face skin
150-11650	Foam set
151-11750	Wig RAAST
150-11850	Cricoid structure
150-13950	Keyboard
151-14050	Chest skin cpl.
150-19550	Technical Manual and block diagram
151-19950	Dir. for Use

**Parts for authorized
service personnel are listed in the
technical manual**



150-10750
Stomach bladder



150-10650
Lung



200-10550
AC Adapter (multi)



282100
Manual defib. adapters
(set of 2)

945004
Flat cable to manikin
(For MRx,
HS 4000, FR and FR2)



321040
Barrel to manikin
(For MRx, HS 4000)



150-13750
Cable assy (USB and DC)




200-11250
Cable, operating device



110100 Carrying case



110104
Wheels, set of 2





200-00550
Cuff assy BP



150-19050
Software CD w/USB cable



310680
Battery holder
(Use 8x 1.5V size D batteries)



150-10750
Stomach bladder



150-10650
Lung



150-10550
Breathing bladder



282100
Manual defib. adapters
(set of 2)



945004
Flat cable to manikin
(For MRx,
HS 4000, FR and FR2)



200-10550
AC Adapter (multi)



153910
Protecting pouch



150-12001 Remote Control FCC-ver.
150-12005 Remote Control (J) FCC-ver.
150-12026 Remote Control (CHN) FCC-ver.
150-12027 Remote Control (Kor) FCC-ver.
150-12033 Remote Control CE-ver.

321040
Barrel to manikin
(For MRx, HS 4000)



150-13750
Cable assy (USB and DC)



150-13850
Cable, BP extension



200-11250
Cable, operating device

NOT ILLUSTRATED

- 150-10450 Air pump
- 150-10950 Clothing
- 150-11250 Upper airway cpl. w/teeth
- 150-11450 Teeth set
- 150-11550 Face skin
- 150-11650 Foam set
- 150-11750 Wig
- 150-11850 Cricoid structure
- 150-13650 Airway closure valve
- 150-13950 Keyboard
- 150-14050 Chest skin cpl.
- 150-14250 Chest skin w/speakers
- 150-19550 Technical Manual and block diagram
- 150-19850 Educational Support
- 150-19950 Dir. for Use
- 312052 Leg assy left w/wounds
- 312053 Leg assy right w/wounds
- 200-00350 Carry case full body, soft

Parts for authorized service personnel are listed in the technical manual



381102
Fastener for Neck Skin
(pkg. 10)

211-10050
Head Assy

211-10150
Head Skin,



00-03050
upil set

381105
Neck Skin (pkg. 6)



211-10250
Airway assy



381106
Teeth Upper Rigid (polyester)



381107
Teeth Upper Soft (vinyl)

200-00250
Cricothyroid Membrane
Tape (1 roll)



380700
IV Arm, Right



380200
Arm Assy,
blood pressure



380210
Cuff Assy, Blood Pressure



380435
Breathing Bladder

211-11050
Leg Assy right
with pedal pulses



300-00150
Thigh pad

211-12050
Leg Assy left
with pedal pulses



380455
Chest skin
with chest drain



312029
Skin & Vein,
IV arm

380461
Kit Adult Female
Genitalia w/valves

380462
Kit Adult Male
Genitalia w/valves



NOT ILLUSTRATED

260305	Power-cord (US)
260306	Power-cord (EUR)
260307	Power-cord (UK)
380100	Link box
380901	Jacket
380902	Trousers
380201	Hardware set, arms to
380400	Chest plate assy
380420	Lung
380446	Foam in chest
380460	Genitalia blank
380470	Belly plate assy
380471	Plug, belly plate
380475	Pin, pelvis
380480	Thrust assy cpl.
211-11650	Cables set legs w/PP
211-11750	Hoses set legs w/PP
211-15150	Foam for valve block,modified
211-15250	SimMan Compression Shaft RP G2005
211-15350	Chest plate assy and shaft G2005
211-19050	Dir. for Use SimMan Mk.2
211-19950	CD-ROM SimMan SW
211-20150	Diabetic foot right w/pedal pulse
211-20250	Varicose vein leg w/pedal pulse
211-25150	Foot right crushed w/pedal pulse
211-25250	Leg assy lower left contused F/A w/PP
211-25350	Leg assy lower right w/fracture w/PP
245-18050	Transportation case Simulator peripherals
245-95050	Monitor, touch screen
245-96050	Web Camera
245-98050	USB Hub
245-99050	PDA
270-00250	IV Bag transfer set with tubing
211-19550	Technical manual SimMan

Advanced Video System

210-31050	Audio/Video adapter (5)
210-31150	SW AVS Video Client & Debrief viewer
210-31250	Audio/Video adapter cable
210-31350	Dir. for Use, AVS

Parts for authorized service personnel are listed in the technical manual



381855
Back pack,
Portability Kit



380810
Pad Set, Shoulder and Hip



383110
Chest Tube Insertion Modules
(pkg. 6)



380410
Post set, ECG/defibrillation



381300
Pulse Oximeter Probe



381201
Jet-Vent Adapter Kit



211-15450 Defib.
studs ext.length



381010
Tubing Assy, Compressor-
Manikin, Air/CO₂



380405
Bladder Assy Mid-Clavicular
(Pneumothorax Chest)



380406
Bladder Set Mid-Axillary
(Pneumothorax Side)



245-16750
Audio cable 1-2



245-16350
Cable Assy, Manikin-Link Box
(15 pin D-sub)



271500
Cable PC/Monitor to Link Box
(9 pin D-sub)

SimMan Cat.No. 380000



381102
Fastener for Neck Skin
(pkg. 10)

381105
Neck Skin (pkg. 6)



381100
Head Assy,
without bronchial
branches



381101
Head Skin,
with airway



381106
Teeth Upper Rigid
(polyester)



381107
Teeth Upper Soft (vinyl)

380110
Remote Control,
Cardiac



380111
Remote Control,
Airway



380700
IV arm,
right



380200
Arm Assy,
blood pressure



380210
Cuff Assy, Blood Pressure

380650
Leg Assy Plain,
right



300-00150
Thigh pad

380600
Leg Assy Plain,
left



380455
Chest skin
with chest drain



312029
Skin & Vein,
IV arm

380461
Kit Adult Female
Genitalia w/valves

380462
Kit Adult Male
Genitalia w/valves



NOT ILLUSTRATED

1005146	Tubing Assy 40ft Compressor-Manikin, Air/CO ₂
1005164	Cable Assy 25 ft, Manikin-Link Box (15-pin D-sub)
260305	Power-cord (US)
260306	Power-cord (EUR)
260307	Power-cord (UK)
380100	Link box
380901	Jacket
380902	Trousers
381900	Directions for Use
380150	Patient Monitor SimMan
380201	Hardware set, arms to
380400	Chest plate assy
380420	Lung
380445	Shaft, compression
380446	Foam in chest
380460	Genitalia blank
380470	Belly plate assy
380471	Plug, belly plate
380475	Pin, pelvis
380480	Thrust assy cpl.
270-00250	IV Bag transfer set with tubing
381950	Technical manual

**Parts for authorized
service personnel are listed in the
technical manual**



380810
Pad Set, Shoulder and Hip



383110
Chest Tube Insertion Modules
(pkg. 6)



381300
Pulse Oximeter Probe



380405
Bladder Assy Mid-Clavicular
(Pneumothorax Chest)



380406
Bladder Set Mid-Axillary
(Pneumothorax Side)



380410
Post set, ECG/defibrilla-
tion



380435
Breathing Bladder



381201
Jet-Vent Adapter Kit



211-15450 Defib.
studs ext.length



381150
Audio Cable



380310
Cable Assy, Manikin-Link Box
(15 pin D-sub)



200-00250
Cricothyroid Membrane
Tape (1 roll)



381010
Tubing Assy,
Compressor-Manikin,
Air/CO₂



271500
Cable PC/Monitor to Link Box
(9 pin D-sub)



200-03050
Pupil set



245-10150
Head skin w/foam



245-13050
Chest skin w/foam



245-12150
Skin and vein IV Arm
right, SimBaby

245-11150
Leg lower IV/IO right
(pkg.3)



245-11650
Leg lower IV/IO left
(pkg.3)



245-19050
Dir. for Use SimBaby



245-13550
Clothing SimBaby

NOT ILLUSTRATED



260305	Power-cord (US)
260306	Power-cord (EUR)
260307	Power-cord (UK)
380100	Link box
245-18050	Transportation case
	Simulator peripherals
245-19550	Technical Manual
245-19950	CD-ROM SimBaby
245-95050	Monitor, touch screen
245-96050	Web Camera
245-98050	USB Hub
245-99050	PDA

Advanced Video System

210-31050	Audio/Video adapter (5)
210-31150	SW AVS Video Client & Debrief viewer
210-31250	Audio/Video adapter cable
210-31350	Dir. for Use, AVS

**Parts for authorized
service personnel are listed in the
technical manual**



220-00850
Face Skin



220-02050
Cyanosis EL panel



220-00450
Baby Hat



220-01050
Right Arm
Skin



220-00950
Torso Skin



220-01150
Left Arm
Skin



220-01950
IO-Fill, Empty
Syringe, 60CC



220-01350
Bladder;
Right Lung



220-01650
Bladder;
Left Lung



220-01250
Bladder; Spontaneous
Breathing



220-01750
Bladder;
Stomach Distension





220-10050
SimNewB
Service DVD



220-29950
SimNewB
CD-ROM Basic



220-19950
SimNewB
Software Advanced



1008531
SimNewB QuickStart Guide

SimNewB™ Preset Patient State Levels

Patient Level	Cyanosis	Heart Rate	Muscle Tone	Respirations	Vocal Sounds	Lung Sounds	Blood Pressure
L5	Off	140/min	Motion	Regular 40/min	Strong Cry	Normal	60/40
L4	On	160/min	Motion	Regular 60/min	Grunting	Coarse Crackles	60/40
L3	On	120/min	Tone	Irregular 10/min	Weak Cry	Coarse Crackles	45/20
L2	On	70/min	Limp	Gasping 4/min	-	Coarse Crackles	30/20
L1	On	40/min	Limp	0	-	-	30/20
L0	On	0	Limp	0	-	-	0/0

220-02250
Level Chart

Cat no 220-02250 rev D




200-10550
AC Adapter, multi



200-11250
Cable operating device



220-03650
Meconium
Assembly





200-01850
Bladder Ass'y, MidAxillary



380405
Bladder Ass'y Mid-Clavicular
(Pneumothorax Chest)



383110
Chest Tube Insertion Modules
(pkg. 6)



380410
Post set, ECG/defibrillation



211-15450 Defib.
studs ext.length



200-00150 Bladder Kit MC Kelly



200-00550 Cuff Ass'y, Blood
Pressure VitalSim



200-01450 Cable-Tube Ass'y,
VitalSim to Manikin



380470 Belly Plate
Ass'y

380471 Plug, Belly Plate




200-03650 Genitalia Kit,
Male & Female



205-01050
Bulb Assembly
with connectors

NOT ILLUSTRATED

200-00350	Case Carrying-Full Body
200-02050	SET HWDR-TOOL-ADLT HIP
200-03750	Chest Foam, MegaCode Kelly VitalSim
200-03850	SET-HDWR-POSTS-TOOLS
205-00150	Shaft, Compression
205-01750	Thigh Tank Hardware
205-10101	Manual Air Pump
300-00250	Pad Ass'y,Ventral Gluteal Injection Male
300-00350	Pad Ass'y, Gluteal Injection Male
300-00450	Valve-Clamp Set, Urinary Hosp
300-04450	Plug,Arm (Left Arm)
380460	Genitalia, Blank
380475	Pin, Pelvis
381102	Strap Set, Head Skin





200-01850
Bladder Assay, Mid-Axillary



380405
Bladder Assay Mid-Clavicular
(Pneumothorax Chest)



383110
Chest Tube Insertion Modules
(pkg. 6)



380410
Post set, ECG/defibrilla-
tion



200-00150 Bladder Kit MC Kelly



211-15450 Defib.
studs ext. length



200-01450 Cable-Tube Assy,
VitalSim to Manikin



380470 Belly Plate
Assy

380471 Plug, Belly Plate




200-03650 Genitalia Kit,
Male & Female



205-01050
Bulb Assembly
with connectors

NOT ILLUSTRATED

200-00350	Case Carrying-Full Body
200-02050	Set hwdr-Tool-adlt Hip
200-00850	Chest Foam, MegaCode Kelly VitalSim
211-15250	Shaft, Compression
300-00250	Pad Assy,Ventral Gluteal Injection Male
300-00350	Pad Assy, Gluteal Injection Male
300-00450	Valve-Clamp Set, Urinary Hosp
380460	Genitalia, Blank
380475	Pin, Pelvis
381102	Strap Set, Head Skin





200-01850
Bladder Assay, Mid-Axillary



380405
Bladder Assay Mid-Clavicular
(Pneumothorax Chest)



383110
Chest Tube Insertion Modules
(pkg. 6)



380410
Post set, ECG/defibrillation



211-15450
Defib. studs ext. length



200-00150 Bladder Kit MC Kelly



200-00550 Cuff Assy, Blood
Pressure VitalSim



200-01450 Cable-Tube Assy,
VitalSim to Manikin



380470 Belly Plate
Assy

380471 Plug, Belly Plate



200-03650 Genitalia Kit,
Male & Female



205-01050
Bulb Assembly
with connectors

NOT ILLUSTRATED

200-00350	Case Carrying-Full Body
200-02050	Set hwdr-Tool-adlt Hip
200-00850	Chest Foam, MegaCode Kelly VitalSim
211-15250	Shaft, Compression
300-00250	Pad Assy,Ventral Gluteal Injection Male
300-00350	Pad Assy, Gluteal Injection Male
300-00450	Valve-Clamp Set, Urinary Hosp
300-04450	Plug,Arm (Left Arm)
380460	Genitalia, Blank
380475	Pin, Pelvis
381102	Strap Set, Head Skin



380410
Post set, ECG/defibrillation



255-00350
Red Bulb Assembly



231-01350 Cable 10 Ft
VitalSim Trainer to
Manikin



231-00250 Head Assy,
MegaCode Kid

375-71001 Arm,
Pediatric Multi-Venous
IV Training Right




231-01250 Plug,
Arm-Pediatric



231-00450 Arm Assy,
Pediatric Left



231-00550 Skin, Torso
MegaCode Kid



375-70150 Kit,
Skin + Vein Assy
Pediatric IV Arm

231-00750 Leg Assy,
Pediatric -Right
W-Intraosseous Leg

231-00101 Lower
Right Leg Pediatric
Intraosseous



231-00850 Leg Assy,
Pediatric-Left



231-00950 Jacket,
Pediatric



231-01050 Pants,
Pediatric



231-01150 Case,
Carry-Pediatric Full
Body

MegaCode Kid Advanced Cat.No. 231-05050



380410
Post set, ECG/defibrillation



255-00350
Red Bulb
Assembly



231-01350 Cable 10 Ft
VitalSim Trainer to
Manikin



231-00250 Head Assy,
MegaCode Kid



375-71001 Arm,
Pediatric Multi-Venous
IV Training Right

231-01250 Plug, Arm-
Pediatric

231-00450 Arm Assy,
Pediatric Left

231-00550 Skin, Torso
MegaCode Kid



375-70150 Kit,
Skin + Vein Assy
Pediatric IV Arm

231-00750 Leg Assy,
Pediatric -Right
W-Intraosseous Leg

231-00101 Lower
Right Leg-Pediatric
Intraosseous

231-00850 Leg Assy,
Pediatric-Left




231-00950 Jacket,
Pediatric



231-01050 Pants,
Pediatric



231-01150 Case,
Carry-Pediatric Full
Body





312029
Skin & Vein,
IV arm



300-01550 Lung Bag Assy (1)



300-00350 Pad Assy,
Gluteal Injection Male



300-00250 Pad
Assy, Ventral Gluteal
Injection Male



380460 Genitalia,
Blank Male



380471 Plug, Belly Plate



300-04050 Pad Assy,
Right Forearm IV



300-04150 Pad Assy,
Left Forearm IV



300-04350 Hardware Set,
Manikin-Pelvis-Torso



300-04250 Hardware Set,
Manikin-Legs-Arms
to Body


NOT ILLUSTRATED

- 300-00450 Valve-Clamp Set,
Anal Hosp
- 300-00550 Valve-Clamp Set,
Urinary Hosp
- 300-03850 Genitalia Assy, Male
- 300-03950 Genitalia Assy, Female
- 300-01850 Teeth, Upper
- 300-01950 Teeth, Lower

255-00350
Red Bulb
Assembly



300-04850 Gown, Hospital-
Back Tie





312029
Skin & Vein,
IV arm



300-00350 Pad Assy,
Gluteal Injection Male



300-00250 Pad
Assy, Ventral Gluteal
Injection Male



380460 Genitalia,
Blank Male



381403 Stoma,
Infected Colostomy



200-00550 Cuff Assy, Blood
Pressure VitalSim



300-01550 Lung Bag Assy (l)



300-04850 Gown, Hospital-
Back Tie



200-01450 Cable-Tube Assy,
VitalSim to Manikin

NOT ILLUSTRATED

- 300-00450 Valve-Clamp Set,
Anal Hosp
- 300-00550 Valve-Clamp Set,
Urinary Hosp
- 300-02150 Belly Plate Assy, Hosp
3 Hole
- 300-03850 Genitalia Assy, Male
- 300-03950 Genitalia Assy, Female
- 300-04950 Stoma, Soft Style
W-Thru Hole
- 380471 Plug, Belly Plate
- 300-01850 Teeth, Upper
- 300-01950 Teeth, Lower



255-00350
Red Bulb
Assembly




300-04250 Hardware Set,
Manikin-Legs-Arms
to Body



300-04350 Hardware Set,
Manikin-Pelvis-Torso







300-04850 Gown, Hospital-Back Tie



300-04250 Hardware Set, Manikin-Legs-Arms to Body

NOT ILLUSTRATED

300-00450	Valve-Clamp Set, Anal Hosp
300-00550	Valve-Clamp Set, Urinary Hosp
325-03250	Genitalia Assy, Male
325-03350	Genitalia Assy, Female
300-01850	Teeth, Upper
300-01950	Teeth, Lower
325-03150	Neck Cradle



380471 Plug, Belly Plate



300-04350 Hardware Set, Manikin-Pelvis-Torso



325-02250 Insert, Breast Foam-Right




325-02150 Insert, Breast Foam-Left

255-00350 Red Bulb Assembly



325-02050 Skin Assy, Chest-Hospital Female



325-00350 Kit, Skin-Vein Adult Female Multi-Venous IV Arm



325-00150 Pad Assy, Ventral Gluteal Injection Female



325-00250 Pad Assy, Gluteal Injection Female



325-02750 Reservoir Assy, Stomach Hospital Female

300-00950 Denture Set, Upper-Lower



325-00850 Head Assy, Female-Advanced Intubation Head



300-00850 Plug, Tracheostomy

325-00950 Lung Assy, Small, 1 each



381402 Pad Assy, Deltoid Injection



300-02250 Reservoir Assy, Colon Hosp

325-01250 Arm Assy, Blood Pressure Right



375-51001 Arm, Female Multi-Venous IV Training Left



300-02150 Belly Plate Assy, Hosp 3 Hole



300-02050 Belly Plate Assy, Hosp w-Speakers



380471 Plug, Belly Plate



381403 Stoma, Infected Colostomy



300-00150 Pad Assy, Thigh Injection



300-04950 Stoma, Soft Style W-Thru Hole

380650 Leg Assy, Right



380600 Leg Assy, Left



380475 Pin, Pelvis



325-01150 Genitalia, Blank-Female



325-02750 Reservoir Assy, Stomach Hospital Female



325-00150 Pad Assy, Ventral Gluteal Injection Female



325-00250 Pad Assy, Gluteal Injection Female



300-04850 Gown, Hospital-Back Tie



300-04250 Hardware Set, Manikin-Legs-Arms to Body



300-04350 Hardware Set, Manikin-Pelvis-Torso

NOT ILLUSTRATED

- 300-00450 Valve-Clamp Set, Anal Hosp
- 300-00550 Valve-Clamp Set, Urinary Hosp
- 325-03250 Genitalia Assy, Male
- 325-03350 Genitalia Assy, Female
- 300-01850 Teeth, Upper
- 300-01950 Teeth, Lower
- 325-03150 Neck Cradle




325-02050 Skin Assy, Chest-Hospital Female



255-00350 Red Bulb Assembly



200-00550 Cuff Assy, Blood Pressure VitalSim



325-00350 Kit, Skin-Vein Female Multi-Venous IV Arm



325-02250 Insert, Breast Foam-Right



325-02150 Insert, Breast Foam-Left



200-01450 Cable-Tube Assy, VitalSim to Manikin

VitalSim Cat.No. 200-100xx



200-10901 Control Unit FCC replacement unit
200-10933 Control Unit CE replacement unit



200-11450 Battery cover operating device

- 200-11001 Operating device (US) FCC
- 200-11005 Operating device (J) FCC
- 200-11026 Operating device (CHN) FCC
- 200-11027 Operating device (KOR) FCC
- 200-11033 Operating device (EUR) CE



200-10150 Carry Case

NOT ILLUSTRATED

Parts for authorized service personnel are listed in the technical manual (200-12133)



200-11350 Battery tray control unit



200-11150 CD Rom VitalSim



200-10250 Cable USB




200-11250 Cable Operating device



200-10550 AC Adapter multi

Nursing Kid Non VitalSim Capable Cat.No. 350-10050



Nursing Kid VitalSim Capable Cat.No. 350-05050



375-70150 Kit,
Skin+Vein Assy
Pediatric IV Arm



255-00350
Red Bulb
Assembly

350-00250 Head Assy,
Pediatric Hospital



350-00650 Pad Assy
Deltoid Injection



350-00350 Arm Assy Left



350-00850 Pad Assy
Gluteal Injection



350-00750 Pad Assy
Thigh Injection



350-00550 Leg Assy Right



350-00450 Leg Assy Left



350-00950 Gown



231-01150 Case,
Carry Pediatric Full
Body



350-00150 Urinary-Anal Valve
350-01050 Pediatric genitalia kit
(Male, Female and
connectors)

NOT ILLUSTRATED



231-01350 Cable 10 Ft
VitalSim Trainer to
Manikin

Nursing Baby Non VitalSim Capable Cat.No. 365-10050


365-02850 Infant Head Complete

365-00301 Infant IV Arm Complete, Right



365-00650 Kit, Sleeve-Vein IV Arm

365-00850 IV Arm Mandrel Assy



365-00201 Infant IO Leg Complete



365-00450 IO Leg Sleeve

365-01250 IO Leg Mandrel Assy



365-02550 Carry Case



365-02450 Tool Kit



365-02250 Lung Bags (2)

365-02350 Band, Lung Retainers (4)



365-03150 Infant Genitalia Kit Complete




240-01250 Baby Pants

365-00101 Infant IV Leg Complete



365-00550 Kit, Sleeve-Vein IV Leg

365-01050 IV Leg Mandrel Assy



NOT ILLUSTRATED

365-01750 Trach plug

Nursing Baby VitalSim Capable Cat.No. 365-05050

365-02850 Infant Head Complete

365-00301 Infant IV Arm Complete, Right



365-00650 Kit, Sleeve-Vein IV Arm

365-00850 IV Arm Mandrel Assy



365-02250 Lung Bags (2)

365-02350 Band, Lung Retainers (4)

365-03150 Infant Genitalia Kit Complete



240-01250 Baby Pants

365-00201 Infant IO Leg Complete



365-00450 IO Leg Sleeve

365-01250 IO Leg Mandrel Assy

365-00101 Infant IV Leg Complete



365-00550 Kit, Sleeve-Vein IV Leg

365-01050 IV Leg Mandrel Assy

365-02550 Carry Case



365-02450 Tool Kit



231-01350 Cable 10 Ft
VitalSim Trainer to
Manikin

NOT ILLUSTRATED

365-01750 Trach plug

Neonatal Resuscitation Baby Cat.No. 240-00001

240-00950
Head Assembly



240-01350
Neonate Chest
Foam Insert



240-00750
Torso Skin



240-00850
Umbilical
Reservoir



240-00350
Right Arm



240-00450
Left Arm



240-00550
Right Leg



240-00650
Left Leg



240-01250
Baby Pants



240-01450
Umbilical Cord
Simulated Blood Set



250-21050
Manikin Lubricant

240-01050
Carry Case





382450
Chest Skin



200-00150
Kit, Bladder
Advanced ALS
Torso



205-01050
Bulb Assembly
with connectors



100-00650
Tank Top, Blue



100-00750
Carry Case



100-00450
Airways (100)



300-00950
Denture Set.
Upper/Lower



255-00350
Red Bulb
Assembly



100-00150
Mouthpieces
(pkg 6)



300-05150
Arm Assy
Right Plain



380800
Arm Assy,
Left Plain



300-02850
Pad Forearm
Left



300-02450
Pad Forearm
Right



300-00250
Pad Assy,Ventral
Gluteal Injection



300-00350
Pad Assy
Gluteal Injection



381402
Pad Assembly,
Deltoid Injection



300-00150
Pad Assy,Thigh
Injection



300-04950
Stoma, Soft Style
W-Thru Hole



380471
Plug, Belly Plate



380650
Leg Assy, Right



380600
Leg Assy Left



300-04850
Gown, Hospital



380460
Genitalia, Blank



300-03950
Adult Female
Genitalia



300-03850
Adult Male
Genitalia



300-00450
Anal Valve,
Clamp Set



300-00550
Urinary Valve,
Clamp Set




300-04350
Hardware set
Manikin-Pelvis-Torso



300-04250
Hardware set
Manikin-Legs-
Arms to Body

NG Tube and Trach Care Trainer Cat.No. 375-00001



300-10150
Head AMH - Int
w/o Lungs



300-00950
Denture Set.
Upper/Lower



300-10250
Chest Plate



255-00350
Red Bulb
Assembly



300-00850
Plug, Trachestomy



300-01550
Lung Bag
Assembly (I)



300-01350
Reservoir Assy,
Stomach Hospital
Adult Male



250-21050
Manikin Lubricant




100-00650
Tank Top, Blue



100-00750
Carry Case





300-02150
Belly Plate Assy
Three Hole with
Stoma



325-02650
Pelvis, No Bridge



380475
Pelvis Pin



300-02250
Adult Colon
Reservoir



325-00150
Adult Female
Ventral Gluteal Pad



325-00250
Adult Female
Dorsal Gluteal Pad



300-00150
Pad Assembly
Thigh Injection



300-03850
Adult Male
Genitalia



300-03950
Adult Female
Genitalia



375-21250
Right Adult
Thigh Stump



375-21150
Left Adult
Thigh Stump



NOT ILLUSTRATED

301-10250 Belly Plate Assy
One Hole

300-00450
Anal Valve,
Clamp Set



300-00550
Urinary Valve,
Clamp Set



280-04050 GUIDE; NEEDLE-TEFLON VIV



280-22250
NEEDLE ASSY
VIRTUAL
PHLEBOTOMY



280-02250
NEEDLE ASSY
VIRTUAL IV



280-03550
STAND;
NEEDLE-VIV



280-03350
ASSY; NEEDLE
STICK-VIV



280-02450
TETHER;
CATHETER-VIV



280-02150
VIRTUAL IV
CATHETER



280-12150
CATHETER
ASSY INFANT
VIRTUAL IV




280-03750
STRAP;
PACKAGING-
HD VIV



280-02650 Power Supply Haptics Dev





NOT ILLUSTRATED

101-10150	Belly Plate, Adult Plain
381402	Pad Assy, Deltoid Injection
276-00650	Pad, Forearm Open Fracture, Right
201-00150	Head Assy, Adult Intubation w-Lung
200-02150	Jacket, blue
200-02250	Pants, blue
200-00350	Case Carry, Full Body

300-01550
Lung Bag
Assembly (I)



276-00550
Belly Plate,
exposed viscera



300-00350
Pad Assy,
Gluteal Injection



300-00250
Pad Assy,
Ventral Gluteal
Injection



276-01150
Chest Plate
Gunshot wound



381503
Pad, Deltoid
Gunshot Large
Caliber



381504
Pad, Deltoid
Gunshot Small
Caliber



381505
Pad, Ventrogluteal
Gunshot Small
Caliber



381506
Pad, Gluteal
Gunshot Large
Caliber



381507
Pad, Thigh
Impaled



381508
Pad, Thigh
Open Fracture



276-00450
Forearm Left,
with burnt hand



276-00350
Forearm Right,
with Lacerated
Dorsum hand



380460
Genitalia, Blank



381502
Foot, Right
Crushed



381510
Right Calf
W-Fracture



381509
Lower Left Leg
Contused Foot-Ankle



261-02001
Mr. Hurt
Head Simulator



250-21050
Manikin Lubricant




300-04250
Hardware set
Manikin-Legs-
Arms to Body



300-04350
Hardware set
Manikin-Pelvis-Torso





100-00150
Mouthpieces
(pkg. 6)



100-00550
Pulse Simulator Assembly



100-00450
Airways (100)



100-10450
Abdominal Thrust Exposed
Viscera



380470 Belly Plate
Assy

380471 Plug, Belly Plate

381503
Pad, Deltoid
Gunshot Large Caliber



381505
Pad, Ventrogluteal
Gunshot Small Caliber



381504
Pad, Deltoid
Gunshot Small Caliber



381506
Pad, Gluteal
Gunshot Large Caliber



381507
Pad, Thigh impaled



381508
Pad, Thigh
Open fracture



300-00750
Blood Simulated,
4oz red



380475
Pelvis Pin



300-04250
Hardware set



276-00350
Forearm
Right W-Lacerated
Dorsum Hand



276-00450
Forearm Left
W-Burnt Hand



381510
Right Calf
w-fracture, Adult



381509
Lower left Leg
contused
foot ankle

NOT ILLUSTRATED

- | | |
|-----------|---|
| 100-10350 | Chest Skin with
Sucking Chest
Wound |
| 200-02150 | Jacket, blue |
| 200-02250 | Pants, blue |
| 200-00350 | Case Carry, Full Body |



101-03650
Deltoid Pad



101-03450
Chest Plate



101-03750
Thigh Pad



101-04250
Coveralls



101-02050
Arm Assy
Right Tuff
Kelly
100+150#



101-02150
Arm Assy
Left Tuff
Kelly
100+150#



101-02250
Arm Assy
Right Tuff
Kelly 175#

101-02350
Arm Assy
Left Tuff
Kelly 175#

101-00650
Leg Assy Right
Tuff Kelly
100+150#



101-00850
Leg Assy Right
Tuff Kelly
175#




101-00750
Leg Assy Left
Tuff Kelly
100+150#

101-00950
Leg Assy Left
Tuff Kelly
175#



101-00150
Waist belt with D-rings



201-00150
Intubation Head
with lungs



100-00550
Pulse Simulator
Assembly



205-01050
Bulb Assembly
with connectors



260-00950
Bladder
Replacement Kit



250-21050
Manikin Lubricant



100-00650
Tank Top



100-00750
Carry Case



Laerdal

Therapy

Products







M3541A
Carry case
for HeartStart MRx,
incl. 3 accessory
pouches/shoulder
strap



948001
Vehicle Wall Mount
HeartStart MRx



948005
Upper fixation bracket



948004
Lower fixation bracket



M3538A
Battery
14.8-volt Lithium Ion



948003
Power cable



948002
DC Power Supply for MRx

94 80 06
Cadex 2 bay charger

94 80 07
Cadex 4 bay charger



40457C 50mm printing paper
(HS 3000QR, HS 4000, MRx) (10 rolls)

40457D 50mm printing paper (HS 3000QR,
HS 4000, MRx) (80 rolls)

989803138171 75mm printing paper for MRx (10rolls)
989803138181 75mm printing paper for MRx (80rolls)



M3725A Test load for use with
M3508A cable (Plug connector)
M1781A Test load for use with
M3507A cable (Barrel connector)
94 50 02 HeartStart FR to Symbio Tester
adapter cable (not for patient use)

Temperature probes for HeartStart MRx (M3536A/M3535A):

- 21078A** Reusable Skin Surface Temperature Probe
- 21082A** Reusable connector cable for disposable temperature probes.
- 21090A** Disposable Oesophageal/Rectal Temperature Probe
1 box = 20 probes. (Requires 21082A)
- 21091A** Disposable Skin Surface Temperature Probe
1 box = 20 probes. (Requires 21082A)



ECG Cables

STANDARD CABLES		ECG 5 LEAD ECG	
	989803147691 Trunk Cable (1,3 meters) (IEC) (standard from Q3/06) <i>The old long trunk cable (2,7m) can be ordered as sparepart (M1949A)</i>		
	M3526A 3-lead ECG patient cable (snap) (1,6 meter) (IEC)	M1971A 5-lead ECG patient cable (grabber) (extremities, + V1) (III til grabber kabel)	
	M3529A 7-lead ECG patient cable (snap) (IEC) (To be used together with M3526A)	M1978A 5-lead ECG patient cable (grabber) – chest (V2 to V6) (To be used together with M1971A)	

SpO₂

M1191A SpO ₂ Adult (>50 kg)		M1194A SpO ₂ Ear clip (>50 kg)	
	M1192A SpO ₂ Child (15-50 kg)		M1195A SpO ₂ Infant (4-15 kg)
	M1941A SpO ₂ Adapter extension cable (2 meter)m	Nelcor Disposable SpO₂ Sensor  M1943A SpO₂ adapter to Nelcor single use probe (M1902/4B) (1,1 meter) M1904B Adult (>50kg) M1903B Pediatric (15-50kg) M1902B infant (4-15kg)	

NBP Antimicrobial Reusable Cuffs (latex free)

			M1598B 1,5 meter NIBP tubing
			M1599B 3 meter NIBP tubing

Reusable Comfort Cuffs (latex free)

M4552B	Easy Care Cuff, Infant (10 - 15 cm Circ.)		M1572A	Pediatric size cuff
M4553B	Easy Care Cuff, Pediatric (14 - 21.5 cm Circ.)		M1573A	Adult small size cuff
M4554B	Easy Care Cuff, Sm Adult (20.5 - 28.5 cm Circ.)		M1574A	Adult size cuff
M4555B	Easy Care Cuff, Adult (27.5 - 36 cm Circ.)		M1575A	Large Adult size cuff
M4556B	Easy Care Cuff, Adult XL (27.5 - 36 cm Circ.)		M1576A	Thigh size cuff
M4557B	Easy Care Cuff, Large Adult (35 - 45 cm Circ.) Cuff cover length 66.5 cm.		M1578A	Multi-Patient Reusable Cuff Kit (M1573A, M1574A, M1575A, M1576A)
M4558B	Easy Care Cuff, Large Adult XL (35 - 45 cm Circ.) Cuff cover length 83 cm			

External Defibrillation Paddles



M3543A
Water resistant external anterior/
anterior paddle set with patient
contact indicator (PCI)

Pads Cable



M3508A
Pads Cable – (Plug connector)



M3507A
Pads Cable – (Barrel connector)



M4763A
Therapy/CPR Cable

Adapters

05-10200
Pads Adapter
(Barrel connector to plug style)

05-10100
Zoll™ Pads Adapter
Can be used with Zoll I200, I400, I600
and M-Series



05-10000
QUIK-COMBO Pads Adapter.
Can be used with Medtronic Physio-
Control LifePak 9, I0C, I1, I2.



Data Management



861276A01
HeartStart Event Review Pro 3.x
(For MRx, FR, FR2, HSI, FRx and HS
4000.)

989803142521
HeartStart 12-lead Transfer Station.
For MRx. Please contact customer
service for language versions.


HeartStart Review Express Connect.
(For MRx, FR, FR2, HSI, FRx and
HS 4000.) Comes in English, French,
German, Italian, Japanese, and Spanish
versions.

89803146981
256MB MRx Data
Card and Tray



Kat. No	Description	CO ₂ Only	O ₂ /CO ₂	Intubated/ LMA or Combitube	Nasal Only	Oral/Nasal	Duration for use	Primary Applications
M2520A	Pediatric Smart CapnoLine™ O ₂		•			•	8H	10-20 kg
M2522A	Adult Smart CapnoLine™ Plus O ₂		•			•	8H	>55 kg
M2524A	Pediatric Smart CapnoLine™	•				•	8H	10-20 kg
M2526A	Adult Smart CapnoLine™ Plus	•				•	8H	>55 kg
M4686A	Adult/ NIV Line -used with full-face or nasal mask CPAP	•			•	•	8H	>545 kg
M1920A	Adult/Pediatric FilterLine® Set	•		•			6H	>2 kg
M1923A	Infant/Neonatal	•		•			72H	>2 kg
M4686A	Adult/NIV Line	•		•			72H	>2 kg

Dual purpose = a probe that enable monitoring of CO₂, and enable supply of O₂.
Humidified = can be used up to 24 hours. All others can be used to 8 hours.



M2520A
Pediatric Smart CapnoLine O₂



M2526A
Adult Smart CapnoLine



M1920A
Adult/Pediatric



M4686A
Adult NIV Line

M2522A
Adult Smart CapnoLine O₂

Therapy (Shock & Pacemaker) electrodes

Electrodes for FR2, MRx, HS 4000



98980315821 I Adult Pads (1 set)

98980315822 I Adult Pads (1 case = 5 sets)

M3713A
Multifunction Electrode Pads
Adult Plus (1 case = 10 sets)



M3717A
Multifunction Electrode Pads
Pediatric Plus (1 case = 5 sets)
(Not for FR2 and FR2+)

Electrodes for MRx, HS 4000 (with barrel style connector)



M3502A Adult Multifunction Pads
(1 case = 10 sets)



M3503A Infant/Child Multifunction
Pads (1 case = 5 sets)

(Not for FR2 and FR2+)

TRAINING AND SIMULATION



905700
Symbio CSI201 (12-lead rhythm
simulator tester) (Barrel connector)

905800
Symbio CSI201 (12-lead rhythm
simulator tester) (Plug connector)

945004
Flat cable to manikin (For MRx,
HS 4000, FR and FR2)



321040
Barrel to manikin
(For MRx, HS 4000)



Q-CPR



M4761A
Q-CPR Compression Sensor



M4762A
Q-CPR Compression Sensor
Adhesive Pads
(10 pack)

HeartStart FR – FR2+ Cat. Nos. M3841A, M3840A

Batteries



M3863A
Long-Life LiMnO₂ Battery. Disposable and recyclable. For FR2 and FR2+ only.



M3848A
FR2+ Rechargeable 12VDC 1.6Ah Li Ion Battery. (This battery is fully charged in 3 hours, to provide 100 shocks (typical) or 5 hours ECG monitoring (typical). This product is compatible with FR2 defibrillators sold after 12/1/01 only. To use with an earlier FR2, your FR2 must be upgraded. Not all older FR2s are eligible.



M3849A
Charger for the Rechargeable Li Ion Battery (M3848A). This product is compatible with FR2 defibrillators sold after 12/1/01 only. To use with an earlier FR2, your FR2 must be upgraded. Not all older FR2s are eligible.



M3874A
ECG Assessment Module
Requires FR2+ with ECG in the display.




989803158211 Adult Pads (1 set)

989803158221 Adult Pads
(1 case = 5 sets)

M3713A
Multifunction Electrode Pads
Adult Plus (1 case = 10 sets)



M3870A Infant/Child
Reduced-Energy
Defibrillator Pads



M3854A
Datacard and tray



941310
Carrycase Semi-Rigid



941311
Carry Case (Hard Pelican)



941309
FR2+ Carry Case Soft

Training Batteries



M3855A
Charger for FR2 Training
& Admin Pack

M3864A
Training & Administration Pack
for FR2 series Defibrillator



07-10900
Training pads



940050xx Laerdal AED Trainer 2



945050
Remote control
for AED Trainer 2



945020
Training battery



945030
Softpack



945090
Training pads,
Laerdal Link
Technology

NOT ILLUSTRATED

- | | |
|--------|---|
| 945021 | Label status |
| 945025 | AED Trainer 2 (pkg. 10) |
| 945026 | Board AED Trainer 2 |
| 945027 | Speaker with cable |
| 945027 | Cable AED Trainer 2 |
| 945055 | Programming kit for AED Trainer 2 (optional software and cable) |
| 945060 | Technical manual |

HeartStart HSI Cat. Nos. M5066A, M5067A

941340
Standard Carry Case
for HeartStart HSI

941341
Slim Carry Case
for HeartStart HSI



NOT ILLUSTRATED

M5074A Infant/child training pads cartridge
M5093A Replacement adult training pads
M5094A Replacement infant/child training pads
M5089A External manikin adapters (10 on a roll)



M5086A HeartStart Trainer

Shock Electrodes



M5071A
Adult SMART Pads
Cartridge (1 Set)



M5072A
Infant/Child SMART
PADS Cartridge (1 set)

942500
Softpack HSI Trainer



Battery

M5070A
(For HSI and FRx) This battery is a disposable, lithium manganese dioxide, long-life battery that will typically last for four years with the defibrillator in standby mode.



M5073A
Adult Training Pads with
Replaceable Seal



HeartStart FRx Cat. Nos. M861304, 861305

Carry Case

941360
FRx carrying case made
of wipeable urethane.



Battery

M5070A
(For HSI and FRx) This battery is a disposable, lithium manganese dioxide, long-life battery that will typically last for four years with the defibrillator in standby mode.



861307A HeartStart FRx Trainer

Shock Electrodes



989803139261
SMART Pads II



989803139311
Infant/child key used in infant/
child defibrillation with the
HeartStart FRx

989803139271
Training Pads II
(1 pair)

989803139291
Replacement
Training Pads TL
(1 pair)



942510
Softpack FRx Trainer

HeartStart HSI and FRx



Data management

PC applications

Basic version free on the web

Download, view and print incident information captured by the AED

HeartStart Review Express Connect

86131IA01

Basic version + save file and email feature. Available in English, French, German, Italian, Spanish and Japanese.

HeartStart Event Review

M3834A Single licence

989803141811 Site wide license

Basic version + save file in database and add Patient and responder data. Configure AED's. Available in English, French, German, Italian, Spanish and Japanese

Event Review Pro

861276 Single licence

861276-A02 3 installs, one CD,

861276-A03 Site-wide license

HeartStart Event Review + 12-lead ECG reports and full Utstein style reporting

HS Configure.

989803143041

Modify configuration settings in the FRx/HSI via IR from PDA

HeartStart FR2, FR2+, HSI and FRx Cat. Nos. M3841A, M3840A



941343

Standby Emergency Wall Cabinet for AED



941314

Bracket without lock



M3857A

Wall Mount Bracket



989803136531

AED Cabinet Basic

HeartStart 4000 Cat. No. M5500B

M5509A
HeartStart 4000
Carrying case



M3502A
Adult multifunction
defib electrodes
(5 pairs)



M3503A
Paediatric multifunction
defib electrodes (5 pairs)



902201
Monitoring electrodes
(10x3 pcs)



M5517A
AC power module



M3507A
Defibrillation pad cable,
Agilent



M1735A
3-lead ECG monitoring
cable



M5516A
Lead acid battery
(pkg. 1)



M5518A
DC power module



40457C
Printing paper
(2 rolls)



M5510A
Data card 8 Mb.
(pkg. 1)



M1191A
SpO₂ transducer

Alternative Defibrillation Pads



989803158211 Adult Pads (1 set)

989803158221 Adult Pads
(1 case = 5 sets)

M3508A
Heartstream pad connector cable
(not illustrated)



M1781A
Defibrillator test load

M5500-89305
User video
PAL



NOT ILLUSTRATED

M5500-91900 User's guide
M3500-90900 Service manual

Parts for authorized service
personnel are listed in the service manual

HeartStart FR *Cat. Nos. 940010, 940020, 940030*



941311
Carry case
(hard, Pelican)



989803158211 Adult Pads (1 set)

989803158221 Adult Pads
(1 case = 5 sets)



BT1
Battery pack

FOR TRAINING



07-10900
Training pads

04-10400
Training card

HeartStart 911 Cat. No. 930000



930200
Accessory pouch

901100
HeartStart (2000/3000 and
911) Battery



902850
In-Unit Battery Charger 110-240VAC



902700
Battery Charger Adapter



902413
HeartStart
Defibrillation
electrodes
(1 pair)



910020
HeartStart Reporter



920620
Patient cable HS 911



920301
Quick Reporter
Charger Paper
(pkg. 5 rolls)



903800
HeartStart Tester

NOT ILLUSTRATED

Parts for authorized
service personnel only

925209	Clock battery
930100	Main board
930110	MCM adapter board
930120	Battery receptacle
930500	Keyboard
930510	Display
930140	Capacitor bank
930125	Patient cable receptacle

For other service parts
see service manual

HeartStart 3000 Cat. No. 920004 - HeartStart 3000QR Cat. No. 920030



920203 (blue)
920212 (red)
Softcase HeartStart 3000



922101 (blue)
922111 (red)
Softcase HeartStart 3000QR



910020
HeartStart Reporter



921100
HeartStart 3000QR
Printer paper
(pkg. 10 rolls)

920301
Quick Reporter
Charger Paper
(pkg. 5 rolls)



920102
Memory Control Module
Plus Manual



903800
HeartStart Tester



924000
Service Module

924110
Check module

NOT ILLUSTRATED

Parts for authorized service personnel only

922000	Upgrade kit HS-3000QR
925005	Battery adapter
925101	Control board
925102	HV board
925103	MCM adapter board
925105	Display
925116	Keyboard HS-3000
925117	Keyboard HS-3000QR
925108	Printer power board
925112	Printer
925208	Patient cable receptacle
925209	Clock battery

For other service parts
see service manual



920602
Patient cable HS 3000
(3-lead shielded)



902413
HeartStart
Defibrillation electrodes
(1 pair)



902420
HeartStart electrodes
(1 pair)



902201
Monitoring
electrodes (10x3 pcs)

Alternative Defibrillation Pads



945005
Patient cable
HeartStart 3000
to Agilent pads



M3502A
Adult multifunction
defib electrodes (5 pairs)

M3503A
Paediatric multifunction
defib electrodes (5 pairs)



945003
Patient cable
HS-3000 to
FR-electrodes



98980315821 | Adult Pads (1 set)

98980315822 | Adult Pads
(1 case = 5 sets)

901100
HeartStart (2000/3000 and
911) Battery



902850
In-Unit Battery Charger | 110-240VAC



902700
Battery Charger Adapter

560200 Patient valve



540105
Disk membranes
(pkg. 10)

540103
Lip valve

851250 Patient valve w/press. reg.



540105
Disk membranes
(pkg. 10)

540103
Lip valve

851252 Pressure
regulator assembly

851350 Patient valve w/press. reg. and lock clip



540105
Disk membranes
(pkg. 10)

540103
Lip valve

851252
Pressure regulator assembly

851103
Lock clips
(pkg. 10)

MASKS



851700
Silicone mask No.2



851600
Silicone mask No. 0/1



851500
Silicone mask No. 00



875200
Multi-Function Mask Cover 4-5



870220
Silicone adult mask 4-5+
w/Multi Function Mask Cover



860220
Silicone child mask 3-4
w/Multi Function Mask Cover



865200
Multi-Function
Mask Cover 3-4

TUBE AND CONNECTORS



871000
Silicone extension
tube (28 cm)



850900
Manometer
connector



850500
Expiration
diverter



875400
Intake/Reservoir
Valve



531901
O₂ reservoir;
2.6 litres



871300
Attachment ring
f/standard head
strap



870400
Head strap w/
attachm. ring



875400
Intake/Reservoir
Valve



551901
O₂ reservoir;
0.6 litre



860420
Compact case cpl., Adult

860410
Compact case cpl., Preterm/
Paediatric



875400 Intake/Reservoir Valve



510103
Cap for LSR Intake
Valve (pkg. 3)



510404
Intake
membranes
(pkg. 10)



871950 Flap valves
(pkg. 2)



51 01 12
O-rings (pkg. 10)

850700
Display case cpl., Preterm

860300
Display case cpl., Paediatric

870600
Display case cpl., Adult




531907
Intake Valve Adapter
(23 mm OD)



870120
Hanging loop



511700
Wall bracket



521100
Wall mount,
Preterm/Paediatric
display case

572000

Wall mount, Adult
display case

NOT ILLUSTRATED


870550 Wall poster
reassembly guide
870950 Directions for use

Laerdal Suction Unit Cat. Nos. 780000, 780010, 780020

REUSABLE CANISTER SYSTEM



784000 Reusable canister




ABBOTT CANISTER SYSTEM



CANISTER HOLDER



CABLES



BEMIS CANISTER SYSTEM (US Only)



BRACKETS

782600 Wall Bracket w/DC-Power Cord

782610 Wall Bracket w/AC Power Cord US

782620 Wall Bracket w/AC-Power Cord EU

782630 Wall Bracket w/AC-Power Cord UK

782640 Wall Bracket w/o Power Cord





780433
Strap for tubing



780434
Hook for tubing strap



780440
External charger kit



780400
Battery



782400
Side pouch



782300
Shoulder strap



782000
Full covering
carrying bag



782010
Semirigid carrying bag


NOT ILLUSTRATED

Parts for authorized service personnel only

780100	Battery cable
780120	Power cable
780130	Handle protection
780140	Selector main rotary switch
780150	Gasket kit cabinet
780160xx	Cabinet cpl.
780170	Vacuum connector
780180	Foam absorbers (6 pcs.)
780190	Battery lid
780240	Power board
780250	MMI board
780270	Battery contact w/spring
780500	LSU pump assy
780510	Piston w/piston gasket
780520	Piston-arm exchange tool
780530	Umbrella valves
780600	LED display
780110	Technical manual
780300	Directions for use
902850	Battery charger (100-240 V AC)
902720	HeartStart Charger adapter (bay for battery)
780900	LSU conversion kit reusable canister
780910	LSU conversion kit Abbott canister

For other service parts
see service manual

781000 Reusable canister



780431
Vacuum canister tubing/
filter tubing
30 cm LSU



781004
Aerosol filter (12 pk)

781005
Aerosol filter (50 pk)



780413
High efficiency
filter kit
Abbott

Laerdal Compact Suction Unit Cat.No. 880020, 880030

884801
1200 ml Reusable
collection canister
(16 each)



884701
800 ml Disposable collection
canister (48 each)

884703
800 ml Disposable collection
canister (6 each)



884905 Connector
tubing (6 pk)



884901 Bacteria filter
(non-sterile) (12 pack)



884904 Elbow



884903
6' Patient tubing



884401
Line cord for US

884402
Hospital grade line
cord (120 VAC)

884403
Line cord for UK

884404
Line cord for
Continental Europe

884405
Line cord for Australia



884400
AC to DC adapter/
charger



884803 Adapter for
Reusable canister
LCSU



884500
12V DC power cord



884600 Carrying case
(softpack w/shoulder
strap)

884601 Shoulder strap
(only)

884702 800ml collection canister package



884701
800 ml Disposable collection
canister (48 each)



884905 Connector
tubing (6 pk)



884903
6' Patient tubing



884901 Bacteria filter
(non-sterile) (12 pack)



884904 Elbow

NOT ILLUSTRATED

**Parts for authorized
service personnel only**

884302	Battery tray
884105	Regulator adjustment spring (only)
884101	Suction manifold/ regulator subassembly Label, LED indicator
884002	Label, LED indicator
884003	Label, suction regulator
884001	Pump housing upper part
884005	Pump housing Lower part
884104	Vacumn adjustment knob
884301	Rechargeable battery
884200	Circuit board for 880020
884201	Circuit board for 880040
884202	High capacity battery conversion kit
884100	Pump motor assembly, complete
884501	Fuse for 12V DC cord, 3A

885100 Directions for use

884802 1200ml collection canister package



884801
1200 ml Reusable collection
canister (16 each)



884905 Connector
tubing (6 pk)



884903
6' Patient tubing






884901 Bacteria filter
(non-sterile) (12 pack)



884904 Elbow

Laerdal Premier Suction Unit *Cat. No. 770000*



NOT ILLUSTRATED

- 770200 Directions for use
- 770210 Service manual



NOT ILLUSTRATED

Parts for authorized service personnel only

- 770101 Gasket kit cabinet
- 770103 Accessories cover
- 770106 Cover battery
- 770107 Internal tube & connector
- 771500 Circuit board cpl.
- 771501 Fuse kit (5A & 800mA)

For other service parts see technical manual

Laerdal Suction Unit *Cat. Nos. 790009, 790010, 790012, 790013, 790017*



NOT ILLUSTRATED

Parts for authorized service personnel only

790503	LSU wing screw
791208	Female charging socket
794001	White lid

For more product information: www.laerdal.com



Laerdal
helping save lives