



Laerdal[®]

helping save lives



Stifneck[®] Select[™] Extrication Collars

Now with
Integrated Hooks
for Nasal
Cannula

Laerdal Spinal Immobilization Solutions

All the adult Stifneck sizes in one easily adjustable, fully assembled collar. Just measure the patient, select the size, adjust, and lock – four simple steps!

The Select collar fits all adult patients and delivers the same superior immobilization you expect from Stifneck collars and, for the smallest patients, the Pedi-Select collar offers an anatomically correct collar for infants and smaller children.



Stifneck[®] Select[™]
Extrication Collar

Stifneck[®] Pedi-Select[™]
Extrication Collar

www.laerdal.com

Stifneck® Select™ and Pedi-Select™ Extrication Collars

Features

Locking Mechanism:

- Four positions: No-Neck, Short, Regular, and Tall
- Clearly marked and easy to identify
- Ensures stability of collar
- Audible response for safety
- Tactile response for safety
- Assures symmetrical adjustment

Rear Panel:

- Molded-in instructions on the back panel
- Unique, large rear panel allows for:
 - c-spine palpation
 - drainage
 - visualization while promoting ventilation for better patient comfort



Tracheal Opening:

- Oversized tracheal opening allows for easy access for advanced airway procedures
- Optimal visualization for assessment
- Easy carotid pulse monitoring

Integrated Cannula Hooks:

- Two positions: left and right
- Properly secures supplemental oxygen devices for fuss-free transportation

Pedi-Select Features:

- Easily identifiable pink strap
- 3 positions: pedi, midi, mini



- 980010 Stifneck Select Extrication Collar
- 980012 Stifneck Select Extrication Collar, Green
- 980020 Stifneck Pedi-Select Extrication Collar

Benefits

- Four adult collars in one/three pediatric collars in one **saves room in bags and ambulance**
- **Easy-to-use sizing measurements** (finger technique)
- Locking mechanism **ensures stability** of collar and symmetrical adjustment
- Directions on collar make training easy and **guarantees instructions always 'on hand'** for rescuer
- In-service training video **minimizes training time and cost**
- Flat storage **requires minimum space**
- Preassembled - assures **speedy application**
- Oversized tracheal opening **allows optimal visualization** for assessment and treatment
- Constructed from **latex-free** quality materials
- **Radiolucent and MRI compatible**
- "Issues in Spinal Care" CDs available upon request



HOSPITAL



EMS



MILITARY

For additional information or to order,
contact your Local Authorized Laerdal Distributor!
www.laerdal.com



Laerdal®
helping save lives