



Laerdal

helping save lives



Providing realistic and affordable training in maternal and newborn care

MamaNatalie

including NeoNatalie



www.laerdal.com/mamanatalie

MamaNatalie birthing simulator

- and her new born baby helps you introduce simulation.

User friendly realism of the simulators and engaging role plays will make the learning sessions efficient and unforgettable, for learners and facilitators alike.

Improving communication skills - Good communication with the mother is essential to assess, interpret and prompt corrective treatment. Strapped to the instructor, MamaNatalie is unique in its ability to facilitate the development of communication skills.

Highly affordable, portable and easy to use - MamaNatalie is the most affordable birthing simulator available. It can be used by any qualified instructor after only a few minutes introduction. It requires no electricity, is highly portable, and quickly made ready for use in any setting. This all makes MamaNatalie particularly suitable for "in situ" simulation training.

Realistic training in control of postpartum haemorrhage - The number one cause of maternal death during childbirth. The instructor can simulate bleeding of up to 1.5 liters, and easily manipulate the condition of the uterus (from boggy to well contracted, and any condition in between) in response to the student's actions.



The simulator is strapped onto the instructor, who manually controls

- Fetal heart sounds
- Delivery of placenta, also incomplete
- Baby positions in birth tract
- Uterine firmness
- Intensity of bleeding and volume of blood loss



MamaNatalie has unique features for communication training and training control of postpartum hemorrhage




Fetal heart rate monitoring



Urine bladder catheterization



Breech delivery



Vacuum delivery




Oxytocin injection



Controlled cord traction



Incomplete placenta





ORDER INFORMATION

MamaNatalie comes with NeoNatalie newborn simulator in either dark or light complexion.
The MamaNatalie Complete model comprises:

- | | |
|--|---|
| <ul style="list-style-type: none"> - MamaNatalie birthing simulator - Placenta w/umbilical cord - Blood concentrate - NeoNatalie Suction - 2 pairs of gloves (1 normal and 1 long) - Fetal stethoscope - Fluid collection tray - Fluid drain - Urine catheter - 20 ml syringe - NeoNatalie skull w/fontanelas - Backpack - Directions for use | <ul style="list-style-type: none"> - NeoNatalie newborn simulator (dark or light complexion) - Squeeze bulbs for simulation of birth cries, spontaneous breathing, palpable umbilical pulse and auscultation of heart sounds - External umbilical cord and 2 umbilical ties - Sheets to simulate towels - Head cap - Storage/carrying pouch - Directions for use |
|--|---|



MamaNatalie comes with a bag for convenient transport

Cat. no.	Title	Includes
340-00233	MamaNatalie Complete (dark)	Includes NeoNatalie newborn simulator with dark complexion
340-00333	MamaNatalie Complete (light)	Includes NeoNatalie newborn simulator with light complexion

Buy One
Gift One

Mira (20) and Shifa survived a severe birth complication thanks to the presence of a well-trained birth attendant



350,000 mothers and nearly 2 million babies die each year during birth. 99 % of these deaths occur in developing countries.

If you buy a MamaNatalie for use in a developed country before end 2012, Laerdal will donate a second unit in your name to a training program in a developing country.